

# BioMidi & BioPlus

## INSTRUCTIONS FOR USE


**MODELS: BioMidi: 425, 625, EF425**

**BioPlus: 500, 600D, 600W, 660D, 660W, 930, 1270, 1400,  
EF600W, EF660W**

Original Instructions for use

Item No.: 765041538  
Revision No.: 20241028  
Language: English

# Quick Guide – BioMidi & BioPlus



## On/Standby

Press the key to turn the cabinet on. Press the key for 6 seconds to switch to standby. The software version of the cabinet will be shown when the cabinet is turned on, followed by the variant and a display test. The cabinet is ready when the temperature is displayed. The cabinet will automatically start a defrost-cycle when turned on, and terminate it again after a system check.

## Setting the temperature

Temperature adjustments are done by holding the key and pressing either or . Confirm the settings by letting go of the keys.

## User menu and alarm settings

| Menu Access  +  →          | ↓   | →          |  |
|----------------------------|-----|------------|--|
|                            | dC* |            | Dry cooling [HO=Off/H1=On]   |
| Local alarm settings       | LAL | LhL [° C]  | Upper alarm limit. Code for activated alarm [A2]                                 |
|                            |     | LLL [° C]  | Lower alarm limit. Code for activated alarm [A3]                                 |
|                            |     | Lhd [min.] | Delay of upper alarm limit   |
|                            |     | LLd [min.] | Delay of lower alarm limit   |
|                            |     | dA On/Off  | Door alarm. Code for activated alarm [A1]. [1=On/0=Off]                          |
|                            |     | dAd [min.] | Delay of door alarm  |
|                            |     | BU On/Off  | Acoustic signal for alarm codes [A1], [A2] and [A3]. [1=On/0=Off]                |
| External alarm settings    | EAL | EhL [° C]  | Upper alarm limit. Code for activated alarm [A4]                                 |
|                            |     | ELL [° C]  | Lower alarm limit. Code for activated alarm [A5]                                 |
|                            |     | Ehd [min.] | Delay of upper alarm limit   |
|                            |     | ELd [min.] | Delay of lower alarm limit   |
|                            |     | dA On/Off  | Door alarm. Code for activated alarm [A1]. [1=On/0=Off]                          |
|                            |     | dAd [min.] | Delay of external door alarm   |
|                            |     | BU On/Off  | Acoustic signal for external alarm codes [A1], [A4], [A5]. [1=On/0=Off]          |
| Offset of sensors          | CAL | cA [° K]   | Calibration of A-sensor. Reference sensor for the refrigeration system           |
|                            |     | cE [° K]   | Calibration of E-sensor. Reference sensor for the display and alarms             |
|                            |     | cF [° K]   | Offset setting for F-sensor. Reference sensor for the low-temperature protection |
| Low-temperature protection | FP  | ACt On/off | Activation/deactivation of low-temperature protection                            |
|                            |     | tES On     | Test of low-temperature protection   |
|                            |     | SEt [° C]  | Setting of the cut-off temperature for the low-temperature protection            |
|                            |     | PrE [...]  | Read-out of the real-time temperature of the F-sensor                            |
|                            | ALL |            | Activation of escorted alarm limits. [FAS]=limits/[ESC]=follows setpoint         |
|                            | dEF |            | Number of defrosts per 24 hours (4 is factory setting)                           |
|                            | dPS |            | Reference sensor for the display (A, E or F)                                     |

## Other Shortcuts

| Keys | Duration    | Function   |
|------|-------------|--|
| +    | > 3 seconds | Start or stop a defrost  |
| +    | > 6 seconds | Activating/deactivating the keypad lock  |
|      | -           | Shows the temperature setpoint value   |
|      | -           | Shows the highest registered temperature spike (since the last reset of the alarm history) |
|      | -           | Shows the lowest registered temperature spike (since the last reset of the alarm history)  |
| +    | > 3 seconds | Clear and reset alarm history  |
| +  + | > 6 seconds | Reset of the set parameters. Restores factory settings                                     |
| +    | > 3 seconds | Access to user menu and alarm settings   |

\* Only ER

### Example: Setting the upper limits for the alarms; LhL

- ↳ Press and hold (P) + (1) until the display shows LAL
- ↳ Press (P) to select LAL. The upper alarm (LhL) limit is now shown in the display
- ↳ Press (P) to select LhL, 25 is now shown in the display
- ↳ Press (-) or (+) to set the desired value
- ↳ Press (P) to return to LAL
- ↳ Press (+) to reach the next level, LLL
- ↳ Lhd, LLd, dA, dAd and BU are located on the same level
- ↳ Leave the user menu by pressing (P) several times until the cabinet's temperature is shown in the display

|             |      |  |
|-------------|------|--|
| Alarm codes | [A1] | Door alarm "dAd" from LAL and/or EAL has been activated        |
|             | [A2] | The upper alarm limit, (LhL) is or has been activated          |
|             | [A3] | The lower alarm limit, (LLL) alarm is or has been activated    |
|             | [A4] | External high alarm EhL is or has been activated (see page 42) |
|             | [A5] | External low alarm ELL is or has been activated (see page 42)  |

### Cancelling an acoustic alarm

Cancelling a door alarm: [A1] Flashes in the display. Press (P) to cancel.

Cancelling a temperature alarm: [A2, A3] Flashes in the display. Press (P) to cancel.

The display will continue to flash if the temperature is outside the alarm limits, and will continue so until the temperature has recovered.

### Read-out of the max./min. temperature

Read the highest recorded temperature inside the cabinet by holding down (+). Read the lowest recorded temperature inside the cabinet by holding down (-).

### Reading the alarm history - Example [A2]

[A2] Flashes in the display. This means that the temperature has exceeded the set value for the upper temperature limit, LhL.

- Press (P) to cancel the [A2]. The display continues to flash, indicating that there is information in the alarm history.
- Press (+), Htt (High temperature time) is shown. Press (P) to see for how long the temperature was above the set alarm limit.
- Press (P) to return to Htt. Press (+) to reach Ht (Highest temperature).
- Press (P) to read the highest recorded temperature during Htt. Press (P) to return to Ht and press (P) again to leave the alarm history function.

The procedure for reading an [A3] alarm is identical, apart from entering the alarm history with the (-) key. When reading out temperatures below set limits, the parameters are Ltt and Lt. A flashing display with no alarm codes indicates that the alarm codes have been canceled, but the alarm system contains information.

### Resetting the alarm temperature and the alarm history

Resetting of the max./min. and the alarm history is done by holding (-) and (+) for more than three seconds. An acoustic signal will be given when reset is complete.

### Sensor read-out and error codes

| Menu Access (P) + (5) →   | ↓   | (P) → [° C]                              | Display code and its message |  |
|---|-----|--|------------------------------|--|
| Sensor for refrigeration system   | P-A | Value on sensor for refrigeration system | F1                           | Error on sensor for refrigeration system |
| Sensor for evaporator   | P-b | Value for evaporator sensor              | F2                           | Error on evaporator sensor               |
| Sensor for condensor  | P-C | Value for condensor sensor               | F3                           | Error on condensor sensor                |
| Sensor for display and alarms   | P-E | Value for display and alarm sensor       | F5                           | Error on sensor for display and alarms   |
| An overheated condensor can be caused by a clogged condensor<br>- Clean the condensor         |     |  | F7                           | Overheated condensor                     |
| Open door indicator.<br>Alarm [A1] will activate if the door is open longer than alarm limits |     |  | -0-                          | Door open                                |

# Table of Content

|  |           |  |            |
|--|-----------|--|------------|
| <b>Quick Guide – BioMidi &amp; BioPlus</b> ..... | <b>2</b>  | <b>Regular maintenance</b> .....             | <b>53</b>  |
| <b>Table of Content</b> .....                    | <b>4</b>  | Cleaning .....                               | 53         |
| <b>Safety</b> .....                              | <b>5</b>  | Door gasket .....                            | 54         |
| Before you proceed .....                         | 5         | <b>General info</b> .....                    | <b>55</b>  |
| <b>Cabinet components</b> .....                  | <b>6</b>  | Responsibility .....                         | 55         |
| BioPlus .....                                    | 6         | Type/number plate .....                      | 56         |
| BioMidi .....                                    | 8         | Defrost water .....                          | 57         |
| <b>Installation</b> .....                        | <b>10</b> | Door self-closing mechanism .....            | 58         |
| Initial setup .....                              | 10        | Door lock .....                              | 59         |
| Anti tilt bracket .....                          | 12        | Access port .....                            | 60         |
| Surroundings .....                               | 13        | Important .....                              | 61         |
| Voltage-free contact .....                       | 15        | Disposal .....                               | 62         |
| Connection to power .....                        | 16        | <b>Datasheet</b> .....                       | <b>63</b>  |
| Equipotential bonding .....                      | 18        | BioMidi 425 .....                            | 63         |
| <b>Options assembly</b> .....                    | <b>20</b> | BioMidi 625 .....                            | 68         |
| Introduction .....                               | 20        | BioPlus 500 .....                            | 73         |
| Shelf carriers .....                             | 21        | BioPlus 600D .....                           | 78         |
| Perforated shelf .....                           | 22        | BioPlus 600W .....                           | 83         |
| Wire shelf .....                                 | 23        | BioPlus 660D .....                           | 91         |
| Carriers for aluminium drawer .....              | 24        | BioPlus 660W .....                           | 96         |
| Aluminium drawer .....                           | 25        | BioPlus EF600W .....                         | 105        |
| Carriers for 40 kg drawer .....                  | 26        | BioPlus EF660W .....                         | 110        |
| 40 kg drawer .....                               | 27        | BioPlus 930 .....                            | 113        |
| How to lock the shelf on a 40 kg drawer .....    | 28        | BioPlus 1270 .....                           | 118        |
| Carriers for 30 kg drawer .....                  | 30        | BioPlus 1400 .....                           | 127        |
| 30 kg drawer .....                               | 31        | <b>Declaration of conformity</b> .....       | <b>136</b> |
| How to lock the 30 kg drawer .....               | 32        | BioMidi .....                                | 136        |
| <b>Start-up</b> .....                            | <b>34</b> | BioMidi – Accessory code 69 .....            | 137        |
| The digital display .....                        | 34        | BioPlus .....                                | 138        |
| Walkthrough of menu .....                        | 35        | BioPlus – Accessory code 69 .....            | 139        |
| Error codes .....                                | 36        | <b>Piping diagram</b> .....                  | <b>140</b> |
| Dry cool .....                                   | 37        | <b>Piping diagram</b> .....                  | <b>140</b> |
| <b>Local alarm settings</b> .....                | <b>38</b> | BioMidi/BioPlus .....                        | 140        |
| Local high alarm .....                           | 38        | BioPlus – With dual refrigeration .....      | 141        |
| Local low alarm .....                            | 38        | <b>Wiring diagram</b> .....                  | <b>142</b> |
| Local high alarm delay .....                     | 39        | BioPlus RF – With SSR .....                  | 142        |
| Local low alarm delay .....                      | 39        | BioPlus ER – With SSR .....                  | 143        |
| On/Off local door alarm .....                    | 40        | BioPlus ER – With SSR .....                  | 143        |
| Delay for local door alarm .....                 | 40        | BioPlus EF600/660 – With SSR .....           | 144        |
| Buzzer – Acoustic local alarms .....             | 41        | BioPlus ER1270/1400 – With SSR .....         | 145        |
| <b>External alarm settings</b> .....             | <b>42</b> | BioPlus RF1270/1400 –                        |            |
| External high alarm .....                        | 42        | With dual compressor – With SSR .....        | 146        |
| External low alarm .....                         | 42        | BioPlus ER1270/1400 –                        |            |
| External high alarm delay .....                  | 43        | With dual compressor – With SSR .....        | 147        |
| External low alarm delay .....                   | 43        | BioPlus RF/EF – With dual compressor         |            |
| On/Off external door alarm .....                 | 44        | – With SSR .....                             | 148        |
| External door alarm delay .....                  | 44        | BioPlus ER1270/1400 – With dual compressor – |            |
| Buzzer – External acoustic settings .....        | 45        | With LTP and SSR .....                       | 149        |
| <b>Parameter settings</b> .....                  | <b>46</b> | BioMidi RR425/625 – With LTP .....           | 150        |
| Sensor offset .....                              | 46        | BioMidi RF425/625 – With SSR .....           | 151        |
| Escorted/set alarm limits .....                  | 48        | BioMidi EF425 – With SSR .....               | 152        |
| Defrosts/24 hours .....                          | 49        | <b>IQ &amp; OQ</b> .....                     | <b>154</b> |
| Display sensor .....                             | 50        | Installation Qualification .....             | 154        |
| Electric low-temperature protection .....        | 51        | Operation Qualification .....                | 154        |
| <b>Ordinary use</b> .....                        | <b>52</b> | <b>PQ</b> .....                              | <b>162</b> |
| Load line .....                                  | 52        | Performance Qualification .....              | 162        |

Copyright © 2006- Gram BioLine, a division of Gram Scientific ApS, Denmark. All rights reserved.

The content of this publication is owned by Gram BioLine, unless otherwise noted, and is protected by Danish and international copyright laws and provisions. Information and images may not be used, copied or transferred without the express written permission of Gram BioLine.

Manufactured by  
**Gram Scientific ApS**  
 Aage Grams Vej 1 · 6500 Vojens · Denmark  
 Tel: +45 73 20 13 00 · Fax: +45 73 20 13 01  
 e-mail: info@gram-bioline.com  
 www.gram-bioline.com





## Before you proceed

*Make sure to read the instructions for use thoroughly before using the cabinet for the first time.  
In the event of need for product support. Do not hesitate to contact us at: [support@gram-bioline.com](mailto:support@gram-bioline.com)*

This instructions for use is intended for the following product series:

### BioMidi and BioPlus

We recommend that you read this instructions for use thoroughly before using the cabinet for the first time. Gram Bioline does not guarantee safe operation if the cabinet is used for anything other than its intended use. Contents of the instructions for use can be subject to change without notice. No part of this instructions for use may be reproduced in any form without expressed written consent of Gram Bioline. Gram Bioline guarantees the cabinet under certain warranty conditions. Gram Bioline is not responsible for any loss or damage of content.

This instructions for use should be considered as an integral part of the cabinet and should be stored close to the cabinet and to be easily accessible. If the instructions for use is lost, please refer to your local distributor or Gram Bioline for a replacement. For current versions of the instructions for use, please go [www.gram-bioline.com](http://www.gram-bioline.com).

## Intended Use

**Gram BioLine BioPlus and BioMidi** refrigerators (RR and ER) and freezers (RF and EF) are designed and manufactured to provide safe and precise conditions for the items stored.

The cabinets are designed for the following operating ranges:

RR: +2/+20 °C

ER: -2/+20 °C

RF: -25/-5 °C

EF (BioPlus): -35/-5 °C

EF (BioMidi): -40/-5 °C

The user must ensure that the cabinet is used in accordance with its intended use. Abnormal use or use conflicting with the intended use or guidelines stipulated in the product documentation can lead to: danger to patient safety, damage to stored items, damage to cabinet and danger to user. Gram BioLine equipment is designed to be used in a system with monitored additional independent alarms to ensure timely reaction to alarms and thereby maximum item safety.

When storing valuable or temperature-sensitive materials or products, it is advisable to employ a continuously monitoring autonomous alarm system. This alarm system should be designed in a manner that allows authorised individuals to promptly detect each alarm state and take the necessary corrective actions.

## Symbols used throughout the instructions for use



**Hazard**



**Risk of fire/flammable materials**



**Risk of electric shock**



**Risk of explosion/explosive materials**



**Risk of material damage**



**Info**



**Risk of personal injury**



**ATEX information**




**Risk of burning/freezing**

# Cabinet components

## BioPlus

*This part describes the main components pertinent to the user.*



### - ATTENTION -

In the event of technical difficulties or malfunction, always contact Gram BioLine technical support or a Gram BioLine authorised service partner.

**1. Shelves and wall rails**

Ensure that the shelves are fixed correctly, before exerting load on them.  
All shelves or drawers must be supported by at least two shelf or drawer supports each.

**2. Door**

Ensure that the door is closed completely after use.  
To minimise fluctuations in temperature, make the door openings as brief as possible.

**3. Door gasket**

Ensure that the gasket is pliable and in good working order.  
Keep the door gasket clean, find instructions in this instructions for use.

**4. Digital display for controller**

Use the display to show the cabinet temperature and, to set the parameters described in this instructions for use.

**5. Base of cabinet**

Ensure cabinets with legs are levelled properly and cabinets with castors are placed on a level surface and locked as mentioned in this instructions for use.

**6. Compressor**

Ensure it is not dented or shows any other signs of damage.

**7. Condenser and fan**

Ensure it is not dented or shows any other signs of damage.

**8. Controller for refrigeration system**

Enclosure for controller, sensors and other parts that monitor and manage the refrigeration system.  
Ensure it is not dented or shows any other signs of damage.

**9. Re-evaporation tray**

Ensure it is not cracked or shows any other signs of damage.  
It is recommended to clean it before applying power to the cabinet for the first time.

**10. Access port**

Used to lead sensors and similar into the cabinet. Ensure that the access port is sealed properly prior to start up.

**11. Terminal box for voltage-free contact**

Used to connect to an external alarm system. Instructions for connection is found in this instructions for use.  
Remember to set the external alarms (EAL).

**12. Defrost water tube**

Outlet for the defrost water coming from the evaporator tray inside the cabinet.  
Ensure it is not damaged or shows signs of damage.

**13. Pressure equalisation valve**

Used to equalise pressure inside the cabinet when opening the door.  
Ensure its not damaged or shows signs of damage.

**14. Equipotential Bonding**

To ensure compliance with ATEX regulations EN 60079-14.  
See installation section for specifications.

**15. Preload cover to access mains terminal connection**

Mains power connection, secured mechanically by the preload cover. Find installation guide in this instructions for use.




**- ATTENTION -**

If parts show signs of damage; do not use the cabinet and contact Gram BioLine or supplier for further aid.

# BioMidi

*This part describes the main components pertinent to the user.*



## - ATTENTION -

In the event of technical difficulties or malfunction, always contact Gram BioLine technical support or a Gram BioLine authorised service partner.

#### 1. Shelves and wall rails

Ensure that the shelves are fixed correctly, before exerting load on them.  
All shelves or drawers must be supported by at least two shelf or drawer supports each.

#### 2. Door

Ensure that the door is closed completely after use.  
To minimise fluctuations in temperature, make the door openings as brief as possible.

#### 3. Door gasket

Ensure that the gasket is pliable and in good working order.  
**Keep the door gasket clean, find instructions in this instructions for use.**

#### 4. Digital display for controller

Use the display to show the cabinet temperature and, to set the parameters described in this instructions for use.

#### 5. Base of cabinet

Ensure cabinets with legs are levelled properly and cabinets with castors are placed on a level surface and locked as mentioned in this instructions for use.

#### 6. Compressor

Ensure it is not dented or shows any other signs of damage.

#### 7. Condenser and fan

Ensure it is not dented or shows any other signs of damage.

#### 8. Controller for refrigeration system

Enclosure for controller, sensors and other parts that monitor and manage the refrigeration system.  
Ensure it is not dented or shows any other signs of damage.

#### 9. Re-evaporation tray

Ensure it is not cracked or shows any other signs of damage.  
It is recommended to clean it before applying power to the cabinet for the first time.

#### 10. Access port

*Used to lead sensors and similar into the cabinet.*  
*Ensure that the access port is sealed properly prior to start up.*

#### 11. Terminal box for voltage-free contact

Used to connect to an external alarm system. Instructions for connection is found in this instructions for use.  
Remember to set the external alarms (EAL).

#### 12. Defrost water tube

Outlet for the defrost water coming from the evaporator tray inside the cabinet.  
Ensure it is not damaged or shows signs of damage.

#### 13. Pressure equalisation valve

Used to equalise pressure inside the cabinet when opening the door.  
Ensure its not damaged or shows signs of damage.

#### 14. Equipotential Bonding

To ensure compliance with ATEX regulations EN 60079-14.  
See installation section for specifications.



#### **- ATTENTION -**

If parts show signs of damage; do not use the cabinet and contact Gram BioLine or supplier for further aid.

# Installation

## Initial setup

*This part of the instructions for use describes how to set up the cabinet.*



- I-1\*:** Due to safety and operating considerations, the cabinet must not be used outdoors.
- I-2\*:** The cabinet should be installed in a dry and sufficiently ventilated area.
- I-3\*:** To ensure efficient operation, the cabinet should not be installed in direct sunlight or close to heat sources.

### **I-4\*:** Ambient operating temperature range

| BioPlus/BioMidi    | Minimum ambient operating temperature | Maximum ambient operating temperature |
|--------------------|---------------------------------------|---------------------------------------|
| RR with solid door | +10 °C                                | +43 °C                                |
| RR with glass door | +10 °C                                | +38 °C                                |
| ER with solid door | +10 °C                                | +43 °C                                |
| ER with glass door | +10 °C                                | +38 °C                                |
| RF with solid door | +10 °C                                | +43 °C                                |
| EF with solid door | +10 °C                                | +30 °C                                |



The cabinet interior must not be exposed to corrosive atmospheres.



Ensure the use of appropriate personal protective equipment such as gloves, when handling the cabinet.



**I-5\*:** Avoid placement of the cabinet in a chloric/acidic environment due to risk of corrosion.



**I-6\*:** The cabinet is shipped with a protective film that should be removed prior to use.




**I-7\*:** Clean the cabinet with a mild soap solution prior to use.




**I-8\*:** The cabinet is only allowed to lay down for very short durations (for instance handling through a doorway). If the cabinet has been laying down, the cabinet must stand up-right for at least 24 hours prior to use. This enables oil in the compressors to flow back into place.

*This part of the instructions for use describes how to adjust legs/castors on the cabinet (BioPlus).*

**I-9\*:** Cabinets equipped with legs should be levelled as shown in the illustration below.



**I-10-11\*:** For cabinets equipped with castors, the floor must be level to ensure stable positioning and safe use. When the cabinet is positioned, the two front casters should be locked.



**- WARNING -**  
**POTENTIAL ELECTROSTATIC CHARGING HAZARD**


Removing protective packaging and film may cause electrostatic discharge. Protective packaging and film shall not be removed in ATEX zones.

## Anti tilt bracket




**I-12\*:** Cabinets with drawers and/or glass door must be secured to a stable vertical surface, ensuring that the cabinet cannot tip over when the drawers are drawn to the outermost position, or the door is open. Brackets for securing is included.

Find the instructions for the tilt bracket below.



**BioPlus**



**BioMidi**




The tilt brackets must be fitted when installing the cabinet, ensuring that the users, surroundings and stored items are not damaged by the cabinet.



## Surroundings

- I-13\*:** The cabinet's back should be placed as close to the wall as possible. Maximum allowed distance between wall and cabinet is 75 mm.



The cabinet is not suited for storing items that emit vapours which, either by themselves or in combination with other chemicals or water, may condensate and corrode the cabinet and its components.



All items in the cabinet that are not encapsulated, or wrapped, should be covered to reduce the risk of corrosion of the cabinet and its components.




### - For Ex environments -

Open containers inside the storage chamber may impact the ATEX zone classification.



**I-14\*:** There must be at least a 30 mm gap between cabinets.



**I-15\*:** Do not cover the upper part of the cabinet if it has a top mounted compressor.



**I-16\*:** Do not use electrical appliances inside the cabinet.

## Voltage-free contact




*This part of the instructions for use covers the voltage-free contact.*

**I-17\*:** The illustration below shows the three connectors for the relay (used ex. in connecting to CTS or other external monitoring systems). The three connections, are respectively. Common, NO and NC.

The moment when voltage is applied the controller draws the relay, this makes it possible for the controller to respond to both high and low alarms, door alarms and power failures. Temperature alarms and door alarms must be configured in the external alarm settings (EAL) before they will activate the voltage-free contact. Find instructions on setting external alarms in the parameter settings section.

The wires that are connected in the connection block for the voltage-free contact, are secured in place by the press-fit plate that is pressed onto the block, thereby also preventing access to the electrical circuit.

Connection of the voltage-free contact should be done by a qualified installer.



**BioMidi**

## Connection to power

*Read the following part thoroughly before connecting the cabinet. Contact a qualified electrician if in doubt.*

### **When setting up in an ordinary scenario that is not subject to regulations for EN 60079-15 zone 2:**

The appliance may be connected in accordance with applicable local heavy current regulations.

### **Please note that there are special regulations for products that are in accordance with EN 60079-15 zone 2 and EN 60079-14: Explosive atmospheres – Electrical installations design, selection and erection.**

The appliance has been manufactured in accordance with EN 60079-15: Electrical apparatus for explosive gas atmospheres – Part 15: Type of protection II 3G Ex nA nC nL IIB Tx Gc . Zone 2 is the applicable zone.

If the appliance is to be installed in a zone 2 environment, specialist personnel should perform the installation, or be consulted beforehand, in order to ensure that the appliance is installed in compliance with the guidelines currently contained in the standard.

- I-19\*:** The cabinet is intended for connection to alternating current. The connection values for voltage (V) and frequency (Hz) are given on the type/number plate.
- I-20-1\*:** The mains terminal is accessed via the preload cover as shown on the top of the control box. Unscrew the preload cover to access the mains supply terminal. Ensure that the preload cover is reinstalled after plugging in the power cord. The leaf spring in the preload cover must engage and preload the plug of the cord as shown in the illustrations on the following page. Ensure that the mains plug is seated completely in the terminal on the cabinet.

The appliance must be connected to the external power supply using a suitable device which mechanically prevents the plug and socket from being separated unintentionally.

- I-21\*:** The connection must be labelled: **"DO NOT SEPARATE WHEN ENERGIZED"**

## – ATTENTION –

Fuses and similar must never be removed or replaced while the appliance is connected to a power source. The electrical terminal box must never be opened while the appliance is connected to a power source.

The compressor's starting equipment must never be dismantled while the appliance is connected to a power source.

Whenever electrical components are dismantled or replaced, the appliance must be moved to an area in which there is no risk of ignition caused by the electrical components or gases contained in the appliance.

Never use the cabinet if the power plug is damaged. The cabinet should be examined by a Gram BioLine service technician in such cases.

When setting up in an ordinary scenario that is not subject to regulations for Zone 2: The appliance may be connected in accordance with applicable local heavy current regulations.

### **In both cases:**

Use a three-wire plug, if the power outlet is intended for a three-wire plug, the lead in green/yellow insulation should be connected to the ground terminal.

Power must be connected via a wall socket. The wall socket should be easily accessible.

All earthing requirements stipulated by the local electricity authorities must be observed. The cabinet plug and wall socket should then have correct earthing. If in doubt, contact your local supplier or authorised electrician.

### **– For Ex environments –**




Special conditions for safe use may apply to this product when installing in an EN 60079-14 environment. Please, see corresponding Ex certificate for specifications.

### **– Technical Support –**

In the event of technical difficulties always contact Gram BioLine technical support or a Gram BioLine authorised service partner. Never dismantle the terminal box or any other electrical component.

*Only applicable for BioPlus.*



# Equipotential bonding

*This part of the instructions for use describes equipotential bonding.*



**I-22-1\*:** For installation in ATEX Cat. 3 Zone 2 areas, it is mandatory to have a equipotential bonding, it is not sufficient to use protective earth through the mains connection.

To secure equipotential bonding of the unit – The mounted external bonding conductor must be used in accordance with national installation requirements e.g. EN 60079-14.


- Mounting of the bonding conductor should be done according to the illustration below.
- Please find location for connection facilities on the back of the cabinet marked with: "Attention – Equipotential bonding".
- The bonding conductor should be at least 4 mm<sup>2</sup> guage thickness.
- Use a ring terminal to ensure adequate bonding.
- Use the supplied M5 machine screw and washer to attach the bonding conductor to the cabinet. Tighten the machine screw to 3.2 Nm.

Bonding of the cabinet is illustrated on this page.





## **– ATTENTION –**

**Please note:** This location is the only manufacturer-approved location for equipotential bonding.





# Options assembly

## Introduction









## Shelf carriers




## Perforated shelf





## Wire shelf




## Carriers for aluminium drawer






## Aluminium drawer





## Carriers for 40 kg drawer




## 40 kg drawer




## How to lock the shelf on a 40 kg drawer




Unlocked shelf



Push down




Push towards cabinet





Locked shelf










## Carriers for 30 kg drawer




## 30 kg drawer




## How to lock the 30 kg drawer



② Unlocked shelf




③ Locked shelf



# Start-up

## The digital display

The digital display depicted below, shows the cabinet's temperature and indicates if the cabinet is connected to a power source. The following chapter explains the conventional operation of the cabinet and alarms for temperature, door etc.



### O-1\*: On/Standby

Press to turn the cabinet on. Press for 6 seconds to switch to standby. The software version of the cabinet will be shown when turning the cabinet on, followed by the software variant and a display test.

The cabinet is ready when the temperature is displayed. The cabinet will automatically start a defrost-cycle when turned on, and terminate it again after a system check.

The cabinet will always commence operation when initially connected to a power supply. For instance after a power outage or when plugging the cabinet in for the first time.

- **Parameter setting**  
Gives access to the cabinet's configurable parameters.
- **Defrost**  
Defrost in progress
- **Dry cool**  
Dry cool in progress (ER-models)
- **Keypad lock**  
Keypad is locked, no access to functions or menus
- **Temperature setting**  
Setting of temperature setpoint and navigation in the menus
- **On/Standby**  
Turn the cabinet on or switch to standby, and navigation in the menus.

### O-2\*: Temperature setting

Temperature adjustments are done by holding and pressing either or . Confirm the settings by letting go of the keys.

#### - ATTENTION -



Make sure the appliance is switched off at the socket before service is performed on electrical parts. It is not sufficient to switch the cabinet to standby on the key, as current will persist in some electrical parts of the cabinet.

#### - WARNING -



DO NOT OPEN, MAINTAIN OR SERVICE IN AN AREA WHEN AN EXPLOSIVE ATMOSPHERE IS PRESENT.

#### - ATTENTION -



High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied.

#### All-round introduction to navigating the menus

Beyond setting the temperature and On/Standby, , , and are used to navigate the menu and set the parameters of the cabinet.





The keys have the following functions in the menu:

- Open a menu step or confirm a set value in the parameter settings.
- Scroll upwards in a given menu or raise a given value in parameter settings (alarm setpoint for instance).
- Scroll downwards in a given menu or lower a given value in parameter settings.
- Go a step back in the menus.















## Walkthrough of menu

The menu below gives a quick overview of the parameter settings for the cabinet.

### User menu

| Menu Access                |  +  → |  |  |  |
|----------------------------|---|---|---|--|
|                            |   | dC*   |   | Dry cooling [HO=Off/H1=On]   |
| Local alarm settings       | LAL   | LhL   | [° C]   | Upper alarm limit. Code for activated alarm [A2]                                 |
|                            |   | LLL   | [° C]   | Lower alarm limit. Code for activated alarm [A3]                                 |
|                            |   | Lhd   | [min.]  | Delay of upper alarm limit   |
|                            |   | LLd   | [min.]  | Delay of lower alarm limit   |
|                            |   | dA  | On/Off  | Door alarm. Code for activated alarm [A1]. [1=On/0=Off]                          |
|                            |   | dAd   | [min.]  | Delay of door alarm  |
|                            |   | bU  | On/Off  | Acoustic signal for alarm codes [A1], [A2] and [A3]. [1=On/0=Off]                |
| External alarm settings    | EAL   | EhL   | [° C]   | Upper alarm limit. Code for activated alarm [A4]                                 |
|                            |   | ELL   | [° C]   | Lower alarm limit. Code for activated alarm [A5]                                 |
|                            |   | Ehd   | [min.]  | Delay of upper alarm limit   |
|                            |   | ELd   | [min.]  | Delay of lower alarm limit   |
|                            |   | dA  | On/Off  | Door alarm. Code for activated alarm [A1]. [1=On/0=Off]                          |
|                            |   | dAd   | [min.]  | Delay of external door alarm   |
|                            |   | bU  | On/Off  | Acoustic signal for external alarm codes [A1], [A4], [A5]. [1=On/0=Off]          |
| Offset of sensors          | CAL   | cA  | [° K]   | Calibration of A-sensor. Reference sensor for the refrigeration system           |
|                            |   | cE  | [° K]   | Calibration of E-sensor. Reference sensor for the display and alarms             |
|                            |   | cF  | [° K]   | Offset setting for F-sensor. Reference sensor for the low-temperature protection |
| Low-temperature protection | FP  | ACt   | On/off  | Activation/deactivation of low-temperature protection                            |
|                            |   | tES   | On  | Test of low-temperature protection   |
|                            |   | SEt   | [° C]   | Setting of the cut-off temperature for the low-temperature protection            |
|                            |   | PrE   | [...]   | Read-out of the real-time temperature of the F-sensor                            |
|                            | ALL   |   |   | Activation of escorted alarm limits. [FAS]=limits/[ESC]=follows setpoint         |
|                            | dEF   |   |   | Number of defrosts per 24 hours (4 is factory setting)                           |
|                            | dPS   |   |   | Reference sensor for the display (A, E or F)                                     |

### Other shortcuts

| Keys  | Duration    | Function   |
|---|-------------|--|
|  +    | > 3 seconds | Start or stop a defrost  |
|  +    | > 6 seconds | Activating/deactivating the keypad lock  |
|    | –           | Shows the temperature setpoint value   |
|    | –           | Shows the highest registered temperature spike (since the last reset of the alarm history) |
|    | –           | Shows the lowest registered temperature spike (since the last reset of the alarm history)  |
|  +    | > 3 seconds | Reset of the alarm history   |
|  +  +  | > 6 seconds | Reset of set parameters. Restores factory settings   |
|  +    | > 3 seconds | Access to the user menu and alarm settings   |



\* Only ER

## Error codes

The following table covers the different error codes that might occur.


| Display code | Explanation  |
|--------------|--|
| - 0 -        | Door is open.  |
| [A1]         | Door alarm "dAd" from LAL and/or EAL has been activated.   |
| [A2]         | Local upper alarm LhL is or has been activated.  |
| [A3]         | Local lower alarm LLL is or has been activated.  |
| [A4]         | External upper alarm EhL is or has been activated.   |
| [A5]         | External lower alarm ELL is or has been activated.   |
| F1           | Error on the main cabinet sensor. The refrigeration system will use an emergency program to make the cabinet run. Temperature stability will be affected. Service is required.   |
| F2           | Error on the evaporator sensor. Service is required.   |
| F3           | Error on the condenser sensor. Service is required.  |
| F4           | Error on condensor sensor 2. Service is required.  |
| F5           | Error on sensor for display and alarm. Service is required.  |
| F7           | F7 indicates that the temperature of the condenser is too high. Turn off the cabinet and check that the condenser is not covered by undesirable items, and insure that the condenser (and possibly filter) is clean. Service is required if the problem is not alleviated. |

## Acknowledge an acoustic alarm

- Cancelling a door alarm: [A1] Flashes in the display. Press  to cancel.
- Cancelling a temperature alarm: [A2, A3] Flashes in the display. Press  to cancel.









The display will continue to flash if the temperature is outside the alarm limits, and will continue so until the temperature has recovered.


## Latching alarms: [A2], [A3], [A4], [A5]

Due to the potential implications of alarms, the red alarm triangle icon will turn on along with the corresponding alarm code will flash in the display. The alarm state will remain on until acknowledged by pressing .

## Reading the alarm history – Example [A2]

[A2] Flashes in the display – This means that the temperature has exceeded the set value for the upper temperature limit, LhL.

- Press  to cancel the [A2]. The display continues to flash, indicating that there is information in the alarm history.
- Press , Htt (High temperature time) is shown, press  to see for how long the temperature was above the set alarm limit.
- Press  to return to Htt. Press  to reach Ht (Highest temperature).
- Press  to read the highest recorded temperature during Htt. Press  to return to Ht and press  again to leave the alarm history function.

The procedure for reading an [A3] alarm is identical, apart from entering the alarm history with . When reading out temperatures below set limits, the parameters are Ltt and Lt.

A flashing display with no alarm codes indicates that the alarm codes have been cancelled, but the alarm history contains information.



## Dry cool

*The following part covers activation/deactivation of the dry cool feature (only applicable to ER models).*

### dC – Setting the dry cool feature

- ↳ Press and hold  $\text{P}$  +  $\text{1}$  until dC is shown in the display
- ↳ Press  $\text{P}$  to choose “dC”
- ↳ Press  $\text{+}$  or  $\text{-}$  to choose between [H1= Off] [H0=On]
- ↳ Press  $\text{P}$  to confirm
- ↳ Leave the user menu by pressing  $\text{0}$ , press several times until the cabinet's temperature is shown in the display.

**Please note:** Dry cool feature reduces the relative humidity in the cabinet, but does not control it. Activation of the dry cool feature can cause greater fluctuations in the cabinet's temperature during defrost.

# Local alarm settings

## Local high alarm Local low alarm

*The following part covers the setting of the upper and lower temperature alarm limits.*

### O-3\*: LhL – Setting the upper alarm limit [° C]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (P) to select "LhL". The upper alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the upper alarm limit
- ↳ Press (P) to confirm the set value
  - The upper alarm limit is now set, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display

### O-4\*: LLL – Setting the lower alarm limit [° C]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (+) to proceed to "LLL"
- ↳ Press (P) to select "LLL". The lower alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the lower alarm limit
- ↳ Press (P) to confirm the set value
  - The lower alarm limit is now set, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

## Local high alarm delay

## Local low alarm delay

*The following part covers the setting of the delay for the local upper and lower temperature alarm limits.*

### O-5\*: Lhd – Setting the delay of the local upper alarm limit [min.]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (+) several times until "Lhd" is shown in the display
- ↳ Press (P) to select "Lhd". The delay of the upper alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the delay of the upper alarm limit
- ↳ Press (P) to confirm the set value
  - The delay of the upper alarm limit is now set, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display

### O-6\*: LLd – Setting the delay of the local lower alarm limit [min.]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (+) several times until "LLd" is shown in the display
- ↳ Press (P) to select "LLd". The delay of the lower alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the delay of the lower alarm limit
- ↳ Press (P) to confirm the set value
  - The delay of the lower alarm limit is now set, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

## On/Off local door alarm

### Delay for local door alarm

*The following part covers the setting of the door alarm and the delay of the door alarm.*

#### O-7\*: dA – Activate/deactivate of local door alarm

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (+) several times until "dA" is shown in the display
- ↳ Press (P) to select "dA"
- ↳ Press (+) or (–) to activate/deactivate the local door alarm [1 = activated/0 = deactivated]
- ↳ Press (P) to confirm the set value
  - The local door alarm is now configured, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display

#### O-8\*: dAd – Setting the delay of the local door alarm [min.]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (+) several times until "dAd" is shown in the display
- ↳ Press (P) to select "dAd". The delay of the local door alarm is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the delay of the local door alarm
- ↳ Press (P) to confirm the set value
  - The delay of the local door alarm is now configured, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

## Buzzer – Acoustic local alarms

*The following part covers the setting of the acoustic local alarms.*

### O-9\*: BU – Activation/deactivation of the acoustic local alarms

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "LAL"
- ↳ Press (P) to select "LAL". "LhL" is now shown in the display
- ↳ Press (+) several times until "BU" is shown in the display
- ↳ Press (P) to select "BU"
- ↳ Press (+) or (–) to activate/deactivate the local acoustic alarms [1 = activated/0 = deactivated]
- ↳ Press (P) to confirm the set value
  - The local acoustic alarms is configured, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

# External alarm settings

## External high alarm External low alarm

*The following part covers the setting of the upper and lower external temperature alarm limits.*

### O-10\*: EhL – Setting the external upper alarm limit [° C]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) several times until "EAL" is shown in the display
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (P) to select "EhL". The external upper alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the external upper alarm limit
- ↳ Press (P) to confirm the set value
  - The external upper alarm limit is now set, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display

### O-11\*: ELL – Setting the external lower alarm limit [° C]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) several times until "EAL" is shown in the display
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (+) to proceed to "ELL"
- ↳ Press (P) to select "ELL". The external lower alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the external lower alarm limit
- ↳ Press (P) to confirm the set value
  - The external lower alarm limit is now set, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

## External high alarm delay

## External low alarm delay

The following parts covers the setting of the delay of the external upper and lower alarms.

### O-12\*: Ehd – Setting the delay of the external upper alarm limit [min.]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "EAL"
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (+) several times until "Ehd" is shown in the display
- ↳ Press (P) to select "Ehd". The external delay of the upper alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the external delay of the upper alarm limit
- ↳ Press (P) to confirm the set value
  - The delay of the external upper alarm limit is now set, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display

### O-13\*: ELd – Setting the delay of the external lower alarm limit [min.]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "EAL"
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (+) several times until "ELd" is shown in the display
- ↳ Press (P) to select "ELd". The delay of the external lower alarm limit is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the delay of the lower alarm limit
- ↳ Press (P) to confirm the set value
  - The delay of the external lower alarm limit is now set, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

## On/Off external door alarm

### External door alarm delay

*The following parts covers the setting and delay, of the external door alarm.*

#### O-14\*: dA – Activation/deactivation of external door alarm

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "EAL"
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (+) several times until "dA" is shown in the display
- ↳ Press (P) to select "dA"
- ↳ Press (+) or (–) to activate/deactivate the external door alarm [1 = activated/0 = deactivated]
- ↳ Press (P) to confirm the set value
  - The external door alarm is now configured, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display

#### O-15\*: dAd – Setting the delay of the external door alarm [min.]

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "EAL"
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (+) several times until "dAd" is shown in the display
- ↳ Press (P) to select "dAd". The delay of the external door alarm is now shown in the display
- ↳ Press (+) or (–) to set the desired value for the delay of the external door alarm
- ↳ Press (P) to confirm the set value
  - The delay of the external door alarm is now configured, proceed to other parameters by pressing (⏻), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏻) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.



## Buzzer – External acoustic settings

*The following part covers the setting of the acoustic external alarms.*

### O-16\*: BU – Activation/deactivation of the acoustic external alarms

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) to proceed to "EAL"
- ↳ Press (P) to select "EAL". "EhL" is now shown in the display
- ↳ Press (+) several times until "BU" is shown in the display
- ↳ Press (P) to select "BU"
- ↳ Press (+) or (–) to activate/deactivate the external acoustic alarms [1 = activated/0 = deactivated]
- ↳ Press (P) to confirm the set value
  - The external acoustic alarms is configured, proceed to other parameters by pressing (⏏), then navigate by using (+) or (–)
- ↳ Leave the user menu by pressing (⏏) several times until the cabinet's temperature is shown in the display



#### – ATTENTION –

High and low temperature alarms set up in the cabinet's controller (including EAL alarms) must be accompanied by additional redundant independent external alarms to ensure maximum item safety.

# Parameter settings

## Sensor offset

The temperature sensors connected to the controller can be offset independently of each other in the parameter settings cAL.

Offset is used in cases where there are deviations in the cabinet's actual operation compared to the display and/or control measurements by independent temperature monitoring.

The cabinet is equipped with an A-sensor and an E-sensor.

**The A-sensor** is used to manage the cabinet's refrigeration system and is fixated in a given position in the cabinet, not in storage space. The location of the A-sensor must not be altered.

**The A-sensor** should be offset if the actual temperature in the cabinet does not match the setpoint, despite taking the hysteresis into consideration. Offset of A-sensor is named "cA".

**The E-sensor** is placed in the cabinet's storage space and can be moved around in the cabinet to get the desired reference point for temperature. The E-sensor is the default display sensor and reference for the alarms. The E-sensor has no effect on control of the refrigeration system.

**The E-sensor** should be offset if the actual temperature in the cabinet's display, provided that the display sensor for reference is the E-sensor, does not match the independent temperature monitoring used for control. Offset of E-sensor is named "cE".

**The F-sensor** is placed inside the storage space, close to the airflow of the cold air exiting the air distribution system. The location of the F-sensor must not be altered as this will have an effect on when the low-temperature protection activates.

**The F-sensor** should be offset if the cut-off temperature for the low-temperature protection, does not match the setpoint temperature for the low-temperature protection. Offset of the F-sensor is named "cF".

### Practical example of offset:

#### Example 1

The temperature in the cabinet is operating colder than the actual setpoint.

With a setpoint of +4 °C, the actual temperature inside the cabinet is between +2 and +4 °C. The desired temperature range is between +3 and +5 °C. This means that "cA", in this case, should be -1.0K, so that the refrigeration system stops 1.0K before and starts 1.0K later than the setpoint normally otherwise would dictate.

#### Example 2

The temperature in the cabinet is operating warmer than the actual setpoint.

With a setpoint of +4 °C, the actual temperature inside the cabinet is between +4 and +6 °C. The desired temperature range is between +3 and +5 °C. This means that "cA", in this case, should be +1.0K, so that the refrigeration system stops 1.0K later and starts 1.0K earlier than the setpoint normally otherwise would dictate.

## Offset of the A-sensor

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "cAL" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "cAL". "cA" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "cA"
- ↳ Press  $\textcircled{+}$  or  $\textcircled{-}$  to offset the A-sensor
- ↳ Press  $\textcircled{P}$  to confirm the set value
  - The A-sensor is now offset, proceed to other parameters by pressing  $\textcircled{U}$ , then navigate by using  $\textcircled{+}$  or  $\textcircled{-}$
- ↳ Leave the user menu by pressing  $\textcircled{U}$  several times the cabinet's temperature is shown in the display

## Offset of the E-sensor

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "cAL" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "cAL". "cA" is shown in the display
- ↳ Press  $\textcircled{+}$  until "cE" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "cE"
- ↳ Press  $\textcircled{+}$  or  $\textcircled{-}$  to offset the E-sensor
- ↳ Press  $\textcircled{P}$  to confirm the set value
  - The E-sensor is now offset, proceed to other parameters by pressing  $\textcircled{U}$ , then navigate by using  $\textcircled{+}$  or  $\textcircled{-}$
- ↳ Leave the user menu by pressing  $\textcircled{U}$  several times until the cabinet's temperature is shown in the display

## Offset of the F-sensor

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "cAL" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "cAL". "cA" is shown in the display
- ↳ Press  $\textcircled{+}$  until "cF" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "cF"
- ↳ Press  $\textcircled{+}$  or  $\textcircled{-}$  to offset the F-sensor
- ↳ Press  $\textcircled{P}$  to confirm the set value
  - The F-sensor is now offset, proceed to other parameters by pressing  $\textcircled{U}$ , then navigate by using  $\textcircled{+}$  or  $\textcircled{-}$
- ↳ Leave the user menu by pressing  $\textcircled{U}$  several times until the cabinet's temperature is shown in the display

---

## Escorted/set alarm limits

*The following part covers the setting of the escorted or set alarm limits.*

### ALL – Setting of escorted/set alarm limits

- ↳ Press and hold  $\textcircled{P}$  +  $\textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until “ALL” is shown in the display
- ↳ Press  $\textcircled{P}$  to select “ALL”
- ↳ Press  $\textcircled{+}$  or  $\textcircled{-}$  to select set or escorted alarm limits
- ↳ Press  $\textcircled{P}$  to confirm the set value
- ↳ Leave the user menu by pressing  $\textcircled{U}$  several times until the cabinet's temperature is shown in the display

“**Set alarm**” is fixed, operating independently from the setpoint. The alarm limits will not change and the selected values will remain regardless of the setpoint being altered.

“**Escorted alarm**” is fixed and locked to the setpoint. The alarm limits will change according to the altered setpoint.

## Defrosts/24 hours

*The following part covers the setting of the amount of defrosts per 24 hours.*

### O-17\*: dEF – Number of defrosts

- ↳ Press and hold **(P)** + **(1)** for more than 3 seconds
- ↳ Press **(+)** several times until "dEF" is shown in the display
- ↳ Press **(P)** to select "dEF"
- ↳ Press **(+)** or **(-)** to set the desired amount of defrosts per 24 hours (factory setting is 4)
- ↳ Press **(P)** to confirm the set value
- ↳ Leave the user menu by pressing **(U)** several times until the cabinet's temperature is shown in the display

**Please note:** It is very important that defrosts should not be set to 0 for a prolonged period of time, as this will reduce the cooling capacity of the cabinet.

## Display sensor

*The following part covers the setting of which the sensor is to be shown in the display*

### O-18\*: dPS – Selection of reference sensor for the display

- ↳ Press and hold (P) + (1) for more than 3 seconds
- ↳ Press (+) several times until "dPS" is shown in the display
- ↳ Press (P) to select "dPS"
- ↳ Press (+) or (-) to select either the A- or E-sensor
- ↳ Press (P) to confirm the set value
- ↳ Leave the user menu by pressing (U) several times until the cabinet's temperature is shown in the display



The dPS only changes the reference sensor for the display, and not the reference sensor for the alarms.



The reference sensor for the refrigeration system is the A-sensor, this cannot be altered.

## Electric low-temperature protection

*The following part covers the electric low-temperature protection*

### FP – Activation/deactivation of low-temperature protection

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "FP" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "FP". "Act" is now shown in the display
- ↳ Press  $\textcircled{P}$  to select "Act"
- ↳ Press  $\textcircled{-}$  or  $\textcircled{+}$  to activate/deactivate [1 = activated/0 = deactivated]
- ↳ Press  $\textcircled{P}$  to confirm the set value
- ↳ Leave the user menu by pressing  $\textcircled{\text{U}}$  several times until the cabinet's temperature is shown in the display

### FP – Setpoint of low-temperature protection

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "FP" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "FP". "Act" is now shown in the display
- ↳ Press  $\textcircled{+}$  several times until "SEt" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "SEt"
- ↳ Press  $\textcircled{-}$  or  $\textcircled{+}$  to select the setpoint temperature for the low-temperature protection
- ↳ Press  $\textcircled{P}$  to confirm the set value
- ↳ Leave the user menu by pressing  $\textcircled{\text{U}}$  several times until the cabinet's temperature is shown in the display

### FP – Test of low-temperature protection

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "FP" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "FP". "Act" is now shown in the display
- ↳ Press  $\textcircled{+}$  to progress to "tES"
- ↳ Press  $\textcircled{P}$  to select "tES" – Test will then be performed
- ↳ Leave the user menu by pressing  $\textcircled{\text{U}}$  several times until the cabinet's temperature is shown in the display

### FP – Temperature of low-temperature protection sensor

- ↳ Press and hold  $\textcircled{P} + \textcircled{1}$  for more than 3 seconds
- ↳ Press  $\textcircled{+}$  several times until "FP" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "FP". "Act" is now shown in the display
- ↳ Press  $\textcircled{+}$  several times until "Pre" is shown in the display
- ↳ Press  $\textcircled{P}$  to select "Pre"
- ↳ Press  $\textcircled{P}$  to show the low-temperature protection sensor temperature
- ↳ Leave the user menu by pressing  $\textcircled{\text{U}}$  several times until the cabinet's temperature is shown in the display

# Ordinary use

## Load line

*The following part shows how items should be placed and stored in the cabinet.*

Keep the marked areas in the cabinets clear of all items (see illustrations below), thereby ensuring adequate air circulation, and therein cooling. Do not place items beneath the lowest shelf bracket.

All products to be stored, that are not wrapped or packed, must be covered in order to avoid unnecessary corrosion of the inner parts of the cabinet.




Items placed on the bottom of the cabinet will cause the air circulation to be impeded, which reduces the cabinets' performance. The illustrations below depict examples of maximum load height of a BioPlus and BioMidi cabinet.




The cabinet interior must not be exposed to corrosive atmospheres.

BioPlus



BioMidi





## Cleaning

*Inadequate cleaning can lead to the cabinet not functioning properly or at all. Icīs. Quōster terios cotīs? Am ere tatuam discern iment, niri iditiur sultor patuam utermantrae nocaessulem inc faucons ulicaelum opopublicut porum et nonia nocastr ibusper isulost imoena ta, niricaediusu moenius.*



The cabinet should be cleaned internally with a mild soap solution (max. 85° C) at suitable intervals and checked thoroughly before it is put into operation again.

Cleaning agents with a pH of  $5 \pm 1$  can be used when a mild soap solution and/or water is used to remove any substance that might damage cabinet components or surfaces. The cleaning agent should be compatible with materials such as steel, alloy, sheet metal, paint, and plastics.


The compressor compartment and in particular the condenser must be kept free from dust and dirt. This is best done with a vacuum cleaner and a brush. The air filters on the condenser and the front panel should be removed and cleaned with warm water (max. 50° C).

It is recommended that the re-evaporation tray is checked regularly for foreign objects and cleaned at least once a year.


Do not flush the compressor compartment and evaporator with water as this may cause short-circuits in the electrical system.

Cleaning agents containing chlorine or compounds of chlorine as well as other corrosive agents, may not be used, as they might cause corrosion to the stainless panels of the cabinet and the evaporator system.

**The location of the condenser for both bottom and top-mounted compressors are illustrated below**



**BioPlus**



**BioMidi**

## Door gasket

*The following part covers the importance of a properly functioning door gasket.*


The door gaskets are an important part of a cabinet. Impaired door gaskets can lead to increased humidity, iced up evaporator (thus reduced cooling capacity), and in some cases, decreased longevity of the cabinet.

It is therefore very important to be aware of the door gasket's condition. Regular inspection is recommended.

The door gasket should be cleaned regularly with a mild soap solution.

If a gasket is to be replaced, please contact your local Gram BioLine distributor.

The illustration below shows the locatin of the door gasket



## Responsibility


*Read the following carefully, for information on technical safety and responsibility on Gram BioLine products.*



### **- WARNING -**

DO NOT OPEN, MAINTAIN OR SERVICE IN AN AREA WHEN AN EXPLOSIVE ATMOSPHERE IS PRESENT.



During servicing make sure the appliance is switched off at the socket before service is performed on the cabinet. It is not sufficient to switch the cabinet to standby on the On/Standby  key, as current will persist in some electrical parts of the cabinet.



Warranty may be void in the event of the cabinet being used for applications other than its intended use, or otherwise not in accordance with the guidelines stipulated in the instructions for use.



Defective parts must be replaced with original parts from Gram BioLine. Gram BioLine can only guarantee functional and safety requirements on the cabinets, if above mentioned is adhered to.



The cabinet should be checked at least once yearly by a Gram BioLine-authorized technician. The refrigeration system and the hermetically sealed compressor require no maintenance. However the condenser requires regular cleaning.





Be aware that cabinets using hydrocarbons (HC) as refrigerant, may require special handling by qualified technicians.

## Type/number plate

If refrigeration fails, first look to see whether the cabinet has been unintentionally switched off, or whether a fuse has blown.

If the cause of failure cannot be found, contact your supplier quoting type and S/N. This information can be found on the type/number plate.



## Defrost water


*The cabinet creates defrost water, that is directed out into a re-evaporation tray at the back of the cabinet.*

Defrost water is led through a tube in the insulation to the re-evaporation tray at the back of the cabinet.



It is recommended that the re-evaporation tray is checked regularly for foreign objects and cleaned accordingly. This shall only be done while the cabinet is turned off.

Be careful not to damage the defrost water tube and the heating element (located in the tray) when cleaning.



---

## Door self-closing mechanism



*Please note: BioMidi and BioPlus cabinets are equipped with doors that have a door self-closing mechanism.*

The door is equipped with a door self-closing mechanism. Open the door up to 90°, and it will shut by itself. Open the door more than 90° and door will remain open.

## Door lock

*The BioMidi and BioPlus cabinets are equipped with a safety lock.*

The illustrations below shows the safety lock.



## Access port

*All BioLine cabinets are equipped with an access port on the back of the cabinets, this can be used to easily fit external sensors and the like.*


The illustration below shows the access port in the BioPlus cabinet. All access ports are constructed in the same fashion, with a conical polystyrene plug (fitted from the back of the cabinet) and a plastic cap (fitted from the inside of the cabinet).



It is very important to refit the polystyrene plug (position 1) and plastic cap (position 2) after mounting a sensor, probe etc. failing to do so can result in lowered performance or malfunction of the cabinet.

Due to the intended use of the BioMidi and BioPlus cabinets, it is critical to ensure a proper seal in the access port.

Access ports are clearly marked "Access port" on the cabinet.





## Important

*In the event of need for product support. Do not hesitate to contact us at: [support@gram-bioline.com](mailto:support@gram-bioline.com)*



## – IMPORTANT –

1. There may be sharp edges on the cabinet housing, compressor room, and interior. Show due diligence when handling the cabinet, negligence of these precautions can lead to injuries.
2. Be cautious about the potential risk of body parts getting wedged in the frame slot between the door and the cabinet when opening or closing the cabinet. Exercise due diligence to avoid any accidents. Failure to take these precautions may result in injuries.
3. Exercise caution to prevent the potential wedging of body parts in the drawer column between the drawers and the interior of the cabinet. Demonstrate due diligence to avoid accidents, as negligence in observing these precautions may lead to injuries.
4. Be particularly vigilant in relation to closing doors with self-close mechanism as these are spring-loaded. Negligence of these precautions can lead to injuries.
5. Unlocked castors can lead to unexpected movements of the cabinet. Lock the castors after installation. Negligence of these precautions can lead to injuries.
6. The re-evaporation tray, re-evaporation tray heating element, pressure pipes and compressor develops considerable heat during operation. Assure yourself that these components are sufficiently tempered before touching. Negligence of these precautions can lead to injuries.
7. The evaporator develops considerable cold during operation. Reassure yourself that the evaporator is sufficiently tempered before touching. Negligence of this precaution may lead to injuries.
8. The fan may cause injury during operation, avoiding touching the fans while the cabinet is connected to the mains. Negligence of these precautions can lead to injuries.
9. No unauthorised modification are allowed.

## Disposal

Electrical and electronic equipment (EEE) contains materials, components and substances that can be dangerous and harmful to human health and the environment if the waste (WEEE) is not disposed of properly.



Contact your local BioLine distributor when the cabinet needs to be disposed of.



Products that are labelled with a “crossed-out wheeled bin” are electric and electronic equipment. The crossed-out wheeled bin symbolises that waste of this type can not be disposed of with unsorted municipal waste, but must be collected separately.



## BioMidi 425

### General data – BioMidi 425

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, alarms, and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm          |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24,5 mm   |
| Gross volume             | 425 litres   |
| Net volume               | 303 litres   |
| Door                     | Left or right hinged   |
| Material interior        | Aluminium/stainless steel, or stainless steel                            |
| Materiale exterior       | White lacquered steel, or stainless steel                                |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                 |
| Dimensions – W x D x H   | 600 x 731 x 1980/2000 mm (RR & RF)<br>622 x 731 x 1980/2000 mm (EF)      |
| Air system               | BioLine ventilated air distribution system                               |
| Defrost system           | Automatic smart defrost with re-evaporation of the defrost water         |
| IP class                 | IP21   |

## BioMidi RR425 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | +2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | K2+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 205 g                       |
| Refrigeration capacity at -10 °C | 279 Watt                    |
| GWP – CO <sub>2</sub> e          | 293.15                      |
| Energy consumption               | 1.19 kWh/24h                |
| Heat emission 100 %              | 217 Watt                    |
| Heat emission default setpoint   | 50 Watt                     |
| Nominal consumption              | 222 Watt                    |
| Start current                    | 10.5A                       |
| Sound level                      | –                           |

## BioMidi RR425 H – with glass door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | +2/+20 °C                   |
| Ambient temperature range        | +10/+38 °C                  |
| Software variant                 | K2+                         |
| K-Value                          | 0.529 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 205 g                       |
| Refrigeration capacity at -10 °C | 279 Watt                    |
| GWP – CO <sub>2</sub> e          | 293.15                      |
| Energy consumption               | 1.87 kWh/24h                |
| Heat emission 100 %              | 218 Watt                    |
| Heat emission default setpoint   | 78 Watt                     |
| Nominal consumption              | 222 Watt                    |
| Start current                    | 10.5A                       |
| Sound level                      | –                           |



## BioMidi RF425 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 200 g                       |
| Refrigeration capacity at -25 °C | 424 Watt                    |
| GWP – CO <sub>2</sub> e          | 784.4                       |
| Energy consumption               | 4.45 kWh/24h                |
| Heat emission 100 %              | 424 Watt                    |
| Heat emission default setpoint   | 185 Watt                    |
| Nominal consumption              | 426 Watt                    |
| Start current                    | 12.6A                       |
| Sound level                      | 44.7 dB(A)                  |

## BioMidi EF425 H - extended freezer – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -40/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 270 g                       |
| Refrigeration capacity at -40 °C | 378 Watt                    |
| GWP – CO <sub>2</sub> e          | 1058.94                     |
| Energy consumption               | 9.40 kWh/24h                |
| Heat emission 100 %              | 549 Watt                    |
| Heat emission default setpoint   | 392 Watt                    |
| Nominal consumption              | 595 Watt                    |
| Start current                    | 23.4A                       |
| Sound level                      | –                           |

## BioMidi RR425 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | +2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | K2+                         |
| K-Value                          | 0.529 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.27 kWh/24h                |
| Heat emission 100 %              | 260 Watt                    |
| Heat emission default setpoint   | 53 Watt                     |
| Nominal consumption              | 258 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 47.2 dB(A)                  |

## BioMidi RR425 G – with glass door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | +2/+20 °C                   |
| Ambient temperature range        | +10/+38 °C                  |
| Software variant                 | K2+                         |
| K-Value                          | 0.47 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.63 kWh/24h                |
| Heat emission 100 %              | 249 Watt                    |
| Heat emission default setpoint   | 68 Watt                     |
| Nominal consumption              | 258 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 47.2 dB(A)                  |



## BioMidi RF425 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 86 g                        |
| Refrigeration capacity at -25 °C | 374 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.2838                      |
| Energy consumption               | 3.61 kWh/24h                |
| Heat emission 100 %              | 341 Watt                    |
| Heat emission default setpoint   | 150 Watt                    |
| Nominal consumption              | 336 Watt                    |
| Start current                    | 13.6A                       |
| Sound level                      | 46.3 dB(A)                  |

## BioMidi EF425 G – extended freezer – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -40/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 105 g                       |
| Refrigeration capacity at -40 °C | 338 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3465                      |
| Energy consumption               | 8.93 kWh/24h                |
| Heat emission 100 %              | 529 Watt                    |
| Heat emission default setpoint   | 372 Watt                    |
| Nominal consumption              | 565 Watt                    |
| Start current                    | 23.5A                       |
| Sound level                      | –                           |

---

## BioMidi 625

---

### General data – BioMidi 625

| Technical specifications | Data  |
|--------------------------|---|
| Connection               | 230 VAC, 50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, alarms and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm         |
| Alarm ports              | Voltage-free contact  |
| Access port              | 1 pc. ø24,5 mm  |
| Gross volume             | 625 litres  |
| Net volume               | 451 litres  |
| Door                     | Left or right hinged  |
| Material interior        | Aluminium/stainless steel or stainless steel                            |
| Materiale exterior       | White lacquered steel or stainless steel                                |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                |
| Dimensions – W x D x H   | 815 x 731 x 1980/2000 mm  |
| Air system               | BioLine ventilated air distribution system                              |
| Defrost system           | Automatic smart defrost with re-evaporation of the defrost water        |
| IP class                 | IP21  |





## BioMidi RR625 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | +2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | K2+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -10 °C | 314 Watt                    |
| GWP – CO <sub>2</sub> e          | 328.9                       |
| Energy consumption               | 1.33 kWh/24h                |
| Heat emission 100 %              | 247 Watt                    |
| Heat emission default setpoint   | 55 Watt                     |
| Nominal consumption              | 244 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |

## BioMidi RR625 H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | +2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | K2+                          |
| K-Value                          | 0.5574 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 230 g                        |
| Refrigeration capacity at -10 °C | 314 Watt                     |
| GWP – CO <sub>2</sub> e          | 328.9                        |
| Energy consumption               | 1.82 kWh/24h                 |
| Heat emission 100 %              | 245 Watt                     |
| Heat emission default setpoint   | 76 Watt                      |
| Nominal consumption              | 244 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |

---

## BioMidi RF625 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 240 g                       |
| Refrigeration capacity at -25 °C | 616 Watt                    |
| GWP – CO <sub>2</sub> e          | 941.28                      |
| Energy consumption               | 4.74 kWh/24h                |
| Heat emission 100 %              | 553 Watt                    |
| Heat emission default setpoint   | 198 Watt                    |
| Nominal consumption              | 558 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | 47.4 dB(A)                  |



## BioMidi RR625 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | +2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | K2+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.30 kWh/24h                |
| Heat emission 100 %              | 255 Watt                    |
| Heat emission default setpoint   | 54 Watt                     |
| Nominal consumption              | 258 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 46.6 dB(A)                  |

## BioMidi RR625 G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | +2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | K2+                          |
| K-Value                          | 0.5574 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T5 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 102 g                        |
| Refrigeration capacity at -10 °C | 389 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.3366                       |
| Energy consumption               | 1.64 kWh/24h                 |
| Heat emission 100 %              | 253 Watt                     |
| Heat emission default setpoint   | 68 Watt                      |
| Nominal consumption              | 258 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | 46.6 dB(A)                   |

---

## BioMidi RF625 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 95 g                        |
| Refrigeration capacity at -25 °C | 512 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3135                      |
| Energy consumption               | 4.03 kWh/24h                |
| Heat emission 100 %              | 445 Watt                    |
| Heat emission default setpoint   | 168 Watt                    |
| Nominal consumption              | 460 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | 48.4 dB(A)                  |

## BioPlus 500

### General data – BioPlus 500

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                        |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24,5 mm   |
| Gross volume             | 500 litres   |
| Net volume               | 365 litres   |
| Door                     | Left or right hinged   |
| Material interior        | Stainless steel  |
| Materiale exterior       | White lacquered steel or stainless steel   |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                               |
| Dimensions – W x D x H   | 600 x 805 x 2025/2275 mm   |
| Air system               | BioLine ventilated air distribution system   |
| Defrost system           | Automatic smart defrost with re-evaporation of the defrost water                       |
| IP class                 | IP21   |

## BioPlus ER500 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 260 g                       |
| Refrigeration capacity at -10 °C | 314 Watt                    |
| GWP – CO <sub>2</sub> e          | 371.8                       |
| Energy consumption               | 1.37 kWh/24h                |
| Heat emission 100 %              | 318 Watt                    |
| Heat emission default setpoint   | 57 Watt                     |
| Nominal consumption              | 304 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |

## BioPlus ER500 H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5199 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 260 g                        |
| Refrigeration capacity at -10 °C | 314 Watt                     |
| GWP – CO <sub>2</sub> e          | 371.8                        |
| Energy consumption               | 1.88 kWh/24h                 |
| Heat emission 100 %              | 324.08 Watt                  |
| Heat emission default setpoint   | 79 Watt                      |
| Nominal consumption              | 304 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |



## BioPlus RF500 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 250 g                       |
| Refrigeration capacity at -25 °C | 616 Watt                    |
| GWP – CO <sub>2</sub> e          | 980.5                       |
| Energy consumption               | 4.78 kWh/24h                |
| Heat emission 100 %              | 519 Watt                    |
| Heat emission default setpoint   | 199 Watt                    |
| Nominal consumption              | 558 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |

## BioPlus ER500 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 95 g                        |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3135                      |
| Energy consumption               | 1.26 kWh/24h                |
| Heat emission 100 %              | 320 Watt                    |
| Heat emission default setpoint   | 53 Watt                     |
| Nominal consumption              | 318 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 46.2 dB(A)                  |

## BioPlus ER500 G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5199 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 95 g                         |
| Refrigeration capacity at -10 °C | 389 W                        |
| GWP – CO <sub>2</sub> e          | 0.3135                       |
| Energy consumption               | 1.67 kWh/24h                 |
| Heat emission 100 %              | 316 Watt                     |
| Heat emission default setpoint   | 70 Watt                      |
| Nominal consumption              | 318 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | 46.2 dB(A)                   |





## BioPlus RF500 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -25 °C | 374 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3036                      |
| Energy consumption               | 3.83 kWh/24h                |
| Heat emission 100 %              | 351 Watt                    |
| Heat emission default setpoint   | 160 Watt                    |
| Nominal consumption              | 336 Watt                    |
| Start current                    | 10.4A                       |
| Sound level                      | 48.7 dB(A)                  |

## BioPlus 600D

### General data – BioPlus 600D

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                        |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24,5 mm   |
| Gross volume             | 600 litres   |
| Net volume               | 432 litres   |
| Door                     | Left or right hinged   |
| Material interior        | Stainless steel  |
| Materiale exterior       | White lacquered steel or stainless steel   |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                               |
| Dimensions – W x D x H   | 695 x 876 x 1875/2125 mm   |
| Air system               | BioLine ventilated air distribution system   |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                           |
| IP class                 | IP21   |



## BioPlus ER600D H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -10 °C | 314 Watt                    |
| GWP – CO <sub>2</sub> e          | 328.9                       |
| Energy consumption               | 1.30 kWh/24h                |
| Heat emission 100 %              | 312 Watt                    |
| Heat emission default setpoint   | 54 Watt                     |
| Nominal consumption              | 306 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |

## BioPlus ER600D H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5189 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 230 g                        |
| Refrigeration capacity at -10 °C | 314 Watt                     |
| GWP – CO <sub>2</sub> e          | 328.9                        |
| Energy consumption               | 1.64 kWh/24h                 |
| Heat emission 100 %              | 309 Watt                     |
| Heat emission default setpoint   | 68 Watt                      |
| Nominal consumption              | 306 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |

---

## BioPlus RF600D H

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -25 °C | 616 Watt                    |
| GWP – CO <sub>2</sub> e          | 902.06                      |
| Energy consumption               | 4.60 kWh/24h                |
| Heat emission 100 %              | 506 Watt                    |
| Heat emission default setpoint   | 192 Watt                    |
| Nominal consumption              | 580 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |



## BioPlus ER600D G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.26 kWh/24h                |
| Heat emission 100 %              | 316 Watt                    |
| Heat emission default setpoint   | 52 Watt                     |
| Nominal consumption              | 320 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 46 dB(A)                    |

## BioPlus ER600D G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5189 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 102 g                        |
| Refrigeration capacity at -10 °C | 389 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.3366                       |
| Energy consumption               | 1.71 kWh/24h                 |
| Heat emission 100 %              | 321 Watt                     |
| Heat emission default setpoint   | 71 Watt                      |
| Nominal consumption              | 320 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | 46 dB(A)                     |

---

## BioPlus RF600D G

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -25 °C | 512 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3036                      |
| Energy consumption               | 4.03 kWh/24h                |
| Heat emission 100 %              | 444 Watt                    |
| Heat emission default setpoint   | 168 Watt                    |
| Nominal consumption              | 482 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | 47.3 dB(A)                  |

## BioPlus 600W

### General data – BioPlus 600W

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                        |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24.5 mm   |
| Gross volume             | 600 litres   |
| Net volume               | 432 litres   |
| Door                     | Left or right hinged   |
| Material interior        | Stainless steel  |
| Materiale exterior       | White lacquered steel or stainless steel   |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                               |
| Dimensions – W x D x H   | 815 x 756 x 1875/2125 mm   |
| Air system               | BioLine ventilated air distribution system   |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                           |
| IP class                 | IP21   |



## BioPlus ER600W H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -10 °C | 314 Watt                    |
| GWP – CO <sub>2</sub> e          | 328.9                       |
| Energy consumption               | 1.52 kWh/24h                |
| Heat emission 100 %              | 315 Watt                    |
| Heat emission default setpoint   | 64 Watt                     |
| Nominal consumption              | 306 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |



## BioPlus ER600W H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5563 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 230 g                        |
| Refrigeration capacity at -10 °C | 314 Watt                     |
| GWP – CO <sub>2</sub> e          | 328.9                        |
| Energy consumption               | 1.89 kWh/24h                 |
| Heat emission 100 %              | 250 Watt                     |
| Heat emission default setpoint   | 79 Watt                      |
| Nominal consumption              | 306 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |





## BioPlus RF600W H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -25 °C | 616 Watt                    |
| GWP – CO <sub>2</sub> e          | 902.06                      |
| Energy consumption               | N/A                         |
| Heat emission 100 %              | N/A                         |
| Heat emission default setpoint   | N/A                         |
| Nominal consumption              | 580 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |

## BioPlus ER600W H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 2 x 215 g                   |
| Refrigeration capacity at -10 °C | 2 x 314 Watt                |
| GWP – CO <sub>2</sub> e          | 614.9                       |
| Energy consumption               | N/A                         |
| Heat emission 100 %              | 65 Watt                     |
| Heat emission default setpoint   | N/A                         |
| Nominal consumption              | 538 Watt                    |
| Start current                    | –                           |
| Sound level                      | –                           |

## BioPlus ER600W H – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5563 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 2 x 215 g                    |
| Refrigeration capacity at -10 °C | 2 x 314 Watt                 |
| GWP – CO <sub>2</sub> e          | 614.9                        |
| Energy consumption               | 2.79 kWh/24h                 |
| Heat emission 100 %              | 499 Watt                     |
| Heat emission default setpoint   | 116 Watt                     |
| Nominal consumption              | 538 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |

## BioPlus RF600W H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 2 x 197 g                   |
| Refrigeration capacity at -25 °C | 2 x 616 Watt                |
| GWP – CO <sub>2</sub> e          | 1545,268                    |
| Energy consumption               | 6.448 kWh/24h               |
| Heat emission 100 %              | 903 Watt                    |
| Heat emission default setpoint   | 269 Watt                    |
| Nominal consumption              | 1072 Watt                   |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |



## BioPlus ER600W G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.35 kWh/24h                |
| Heat emission 100 %              | 322.75 Watt                 |
| Heat emission default setpoint   | 56 Watt                     |
| Nominal consumption              | 320 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 44.8 dB(A)                  |

## BioPlus ER600W G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5563 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 102 g                        |
| Refrigeration capacity at -10 °C | 389 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.3366                       |
| Energy consumption               | 1.65 kWh/24h                 |
| Heat emission 100 %              | 317.18 Watt                  |
| Heat emission default setpoint   | 69 Watt                      |
| Nominal consumption              | 320 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | 44.8 dB(A)                   |

---

## BioPlus RF600W G

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -25 °C | 512 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3036                      |
| Energy consumption               | 4.07 kWh/24h                |
| Heat emission 100 %              | 444 Watt                    |
| Heat emission default setpoint   | 169 Watt                    |
| Nominal consumption              | 482 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | 45.8 dB(A)                  |



## BioPlus ER600W G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 85 g                    |
| Refrigeration capacity at -10 °C | 2 x 389 Watt                |
| GWP – CO <sub>2</sub> e          | 0.627                       |
| Energy consumption               | 1.953 kWh/24h               |
| Heat emission 100 %              | 561 Watt                    |
| Heat emission default setpoint   | 81 Watt                     |
| Nominal consumption              | 572 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | –                           |



## BioPlus ER600W G – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5563 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 2 x 85 g                     |
| Refrigeration capacity at -10 °C | 2 x 389 Watt                 |
| GWP – CO <sub>2</sub> e          | 0.627                        |
| Energy consumption               | 2.586 kWh/24h                |
| Heat emission 100 %              | 550.5 Watt                   |
| Heat emission default setpoint   | 108 Watt                     |
| Nominal consumption              | 572 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | –                            |

---

## BioPlus RF600W G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 100 g                   |
| Refrigeration capacity at -25 °C | 2 x 512 Watt                |
| GWP – CO <sub>2</sub> e          | 0.66                        |
| Energy consumption               | 5.815 kWh/24h               |
| Heat emission 100 %              | 806 Watt                    |
| Heat emission default setpoint   | 242 Watt                    |
| Nominal consumption              | 876 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | –                           |

## BioPlus 660D

### General data – BioPlus 660D

| Technical specifications | Data  |
|--------------------------|---|
| Connection               | 230 VAC/50 Hz   |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and calibration function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                             |
| Alarm ports              | Voltage-free contact  |
| Access port              | 1 pc. ø24,5 mm  |
| Gross volume             | 660 litres  |
| Net volume               | 484 litres  |
| Door                     | Left or right hinged  |
| Material interior        | Stainless steel   |
| Materiale exterior       | White lacquered steel or stainless steel  |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                                    |
| Dimensions – W x D x H   | 695 x 876 x 2025/2275 mm  |
| Air system               | BioLine ventilated air distribution system  |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                                |
| IP class                 | IP21  |

## BioPlus ER660D H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -10 °C | 314 Watt                    |
| GWP – CO <sub>2</sub> e          | 328.9                       |
| Energy consumption               | 1.32 kWh/24h                |
| Heat emission 100 %              | 313 Watt                    |
| Heat emission default setpoint   | 55 Watt                     |
| Nominal consumption              | 306 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |

## BioPlus ER660D H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5223 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 230 g                        |
| Refrigeration capacity at -10 °C | 314 Watt                     |
| GWP – CO <sub>2</sub> e          | 328.9                        |
| Energy consumption               | 1.84 kWh/24h                 |
| Heat emission 100 %              | 310 Watt                     |
| Heat emission default setpoint   | 77 Watt                      |
| Nominal consumption              | 306 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |





## BioPlus RF660D H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -25 °C | 616 Watt                    |
| GWP – CO <sub>2</sub> e          | 902.06                      |
| Energy consumption               | 4.91 kWh/24h                |
| Heat emission 100 %              | 505 Watt                    |
| Heat emission default setpoint   | 205 Watt                    |
| Nominal consumption              | 580 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |

## BioPlus ER660D G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.18 kWh/24h                |
| Heat emission 100 %              | 317 Watt                    |
| Heat emission default setpoint   | 49 Watt                     |
| Nominal consumption              | 320 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 45.9 dB(A)                  |

## BioPlus ER660D G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5223 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 102 g                        |
| Refrigeration capacity at -10 °C | 389 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.3366                       |
| Energy consumption               | 1.69 kWh/24h                 |
| Heat emission 100 %              | 317 Watt                     |
| Heat emission default setpoint   | 70 Watt                      |
| Nominal consumption              | 320 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | 45.9 dB(A)                   |



## BioPlus RF660D G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -25 °C | 512 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3036                      |
| Energy consumption               | 4.23 kWh/24h                |
| Heat emission 100 %              | 444 Watt                    |
| Heat emission default setpoint   | 176 Watt                    |
| Nominal consumption              | 482 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | 45.6 dB(A)                  |

## BioPlus 660W

### General data – BioPlus 660W

| Technical specifications | Data  |
|--------------------------|---|
| Connection               | 230 VAC/50 Hz   |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and calibration function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                             |
| Alarm ports              | Voltage-free contact  |
| Access port              | 1 pc. ø24,5 mm  |
| Gross volume             | 660 litres  |
| Net volume               | 484 litres  |
| Door                     | Left or right hinged  |
| Material interior        | Stainless steel   |
| Materiale exterior       | White lacquered steel or stainless steel  |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                                    |
| Dimensions – W x D x H   | 815 x 756 x 2025/2275 mm  |
| Air system               | BioLine ventilated air distribution system  |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                                |
| IP class                 | IP21  |



## BioPlus ER660W H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -10 °C | 314 Watt                    |
| GWP – CO <sub>2</sub> e          | 328.9                       |
| Energy consumption               | 1.42 kWh/24h                |
| Heat emission 100 %              | 313 Watt                    |
| Heat emission default setpoint   | 59 Watt                     |
| Nominal consumption              | 306 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |

## BioPlus ER660W H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5604 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 230 g                        |
| Refrigeration capacity at -10 °C | 314 Watt                     |
| GWP – CO <sub>2</sub> e          | 328.9                        |
| Energy consumption               | 1.90 kWh/24h                 |
| Heat emission 100 %              | 307 Watt                     |
| Heat emission default setpoint   | 79 Watt                      |
| Nominal consumption              | 306 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |

## BioPlus RF660W H

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -25 °C | 616 Watt                    |
| GWP – CO <sub>2</sub> e          | 902.06                      |
| Energy consumption               | 5.07 kWh/24h                |
| Heat emission 100 %              | 497 Watt                    |
| Heat emission default setpoint   | 211 Watt                    |
| Nominal consumption              | 580 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |



## BioPlus ER660W H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 2 x 215 g                   |
| Refrigeration capacity at -10 °C | 2 x 314 Watt                |
| GWP – CO <sub>2</sub> e          | 614.9                       |
| Energy consumption               | 2.26 kWh/24h                |
| Heat emission 100 %              | 500 Watt                    |
| Heat emission default setpoint   | 94 Watt                     |
| Nominal consumption              | 538 Watt                    |
| Start current                    | 11.4A                       |
| Sound level                      | –                           |



## BioPlus ER660W H – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5604 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 2 x 215 g                    |
| Refrigeration capacity at -10 °C | 2 x 314 Watt                 |
| GWP – CO <sub>2</sub> e          | 614.9                        |
| Energy consumption               | N/A                          |
| Heat emission 100 %              | 65                           |
| Heat emission default setpoint   | N/A                          |
| Nominal consumption              | 538 Watt                     |
| Start current                    | 11.4A                        |
| Sound level                      | –                            |

---

## BioPlus RF660W H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 2 x 197 g                   |
| Refrigeration capacity at -25 °C | 2 x 616 Watt                |
| GWP – CO <sub>2</sub> e          | 1545.268                    |
| Energy consumption               | 6.77 kWh/24h                |
| Heat emission 100 %              | 963 Watt                    |
| Heat emission default setpoint   | 282 Watt                    |
| Nominal consumption              | 1072 Watt                   |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |





## BioPlus ER660W G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 102 g                       |
| Refrigeration capacity at -10 °C | 389 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.3366                      |
| Energy consumption               | 1.28 kWh/24h                |
| Heat emission 100 %              | 316 Watt                    |
| Heat emission default setpoint   | 53 Watt                     |
| Nominal consumption              | 320 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | 44.8 dB(A)                  |

## BioPlus ER660W G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5604 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 102 g                        |
| Refrigeration capacity at -10 °C | 389 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.3366                       |
| Energy consumption               | 1.78 kWh/24h                 |
| Heat emission 100 %              | 320 Watt                     |
| Heat emission default setpoint   | 74 Watt                      |
| Nominal consumption              | 320 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | 44.8 dB(A)                   |

---

## BioPlus RF660W G

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25 °C/-50 °C               |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -25 °C | 512 Watt                    |
| GWP – CO <sub>2</sub> e          | N/A                         |
| Energy consumption               | 4.46 kWh/24h                |
| Heat emission 100 %              | 447 Watt                    |
| Heat emission default setpoint   | 186 Watt                    |
| Nominal consumption              | 482 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | 49.3 dB(A)                  |



## BioPlus ER660W G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 85 g                    |
| Refrigeration capacity at -10 °C | 2 x 389 Watt                |
| GWP – CO <sub>2</sub> e          | 0.627                       |
| Energy consumption               | 2.20 kWh/24h                |
| Heat emission 100 %              | 551 Watt                    |
| Heat emission default setpoint   | 92 Watt                     |
| Nominal consumption              | 572 Watt                    |
| Start current                    | 9.7A                        |
| Sound level                      | –                           |

## BioPlus ER660W G – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5604 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 2 x 85 g                     |
| Refrigeration capacity at -10 °C | 2 x 389 Watt                 |
| GWP – CO <sub>2</sub> e          | 0.627                        |
| Energy consumption               | 2.175 kWh/24h                |
| Heat emission 100 %              | 556 Watt                     |
| Heat emission default setpoint   | 113 Watt                     |
| Nominal consumption              | 572 Watt                     |
| Start current                    | 9.7A                         |
| Sound level                      | –                            |

## BioPlus RF660W G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25 °C/-5 °C                |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 100 g                   |
| Refrigeration capacity at -25 °C | 2 x 512 Watt                |
| GWP – CO <sub>2</sub> e          | 0.66                        |
| Energy consumption               | 6.02 kWh/24h                |
| Heat emission 100 %              | 821 Watt                    |
| Heat emission default setpoint   | 251 Watt                    |
| Nominal consumption              | 876 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | –                           |

## BioPlus EF600W

### General data – BioPlus EF600W

| Technical specifications | Data  |
|--------------------------|---|
| Connection               | 230 VAC/50 Hz   |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm           |
| Alarm ports              | Voltage-free contact  |
| Access port              | 1 pc. ø24,5 mm  |
| Gross volume             | 600 litres  |
| Net volume               | 432 litres  |
| Door                     | Left or right hinged  |
| Material interior        | Stainless steel and inner doors in PMMA                                   |
| Materiale exterior       | White lacquered steel or stainless steel                                  |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                  |
| Dimensions – W x D x H   | 837 x 756 x 1875/2125 mm  |
| Air system               | BioLine ventilated air distribution system                                |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water              |
| IP class                 | IP21  |

## BioPlus EF600W H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 360 g                       |
| Refrigeration capacity at -40 °C | 313 Watt                    |
| GWP – CO <sub>2</sub> e          | 1411.92                     |
| Energy consumption               | 9.13 kWh/24h                |
| Heat emission 100 %              | 477 Watt                    |
| Heat emission default setpoint   | 380 Watt                    |
| Nominal consumption              | 474 Watt                    |
| Start current                    | 20A                         |
| Sound level                      | –                           |

## BioPlus EF600W H – water cooled, with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 230 g                       |
| Refrigeration capacity at -40 °C | 355 Watt                    |
| GWP – CO <sub>2</sub> e          | 902,06                      |
| Energy consumption               | 7.93 kWh/24h                |
| Heat emission 100 %              | 523 Watt                    |
| Heat emission default setpoint   | 330 Watt                    |
| Nominal consumption              | 494 Watt                    |
| Start current                    | 20A                         |
| Sound level                      | –                           |



## BioPlus EF600W H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 2 x 197 g                   |
| Refrigeration capacity at -40 °C | 2 x 616 Watt                |
| GWP – CO <sub>2</sub> e          | 1545,268                    |
| Energy consumption               | 12.70 kWh/24h               |
| Heat emission 100 %              | 830 Watt                    |
| Heat emission default setpoint   | 529 Watt                    |
| Nominal consumption              | 748 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |

## BioPlus EF600W G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/(m2*K)               |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -40 °C | 244 Watt                    |
| GWP – CO2e                       | 0.3036                      |
| Energy consumption               | 8.86 kWh/24h                |
| Heat emission 100 %              | 442 Watt                    |
| Heat emission default setpoint   | 369 Watt                    |
| Nominal consumption              | 402 Watt                    |
| Start current                    | 19.5A                       |
| Sound level                      | –                           |

## BioPlus EF600W G – water cooled, with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/(m2*K)               |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | –                           |
| Refrigerant charge               | N/A                         |
| Refrigeration capacity at -40 °C | –                           |
| GWP – CO2e                       | N/A                         |
| Energy consumption               | N/A                         |
| Heat emission 100 %              | N/A                         |
| Heat emission default setpoint   | N/A                         |
| Nominal consumption              | N/A                         |
| Start current                    | –                           |
| Sound level                      | –                           |





## BioPlus EF600W G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 94 g                    |
| Refrigeration capacity at -40 °C | 2 x 286 Watt                |
| GWP – CO <sub>2</sub> e          | 0.594                       |
| Energy consumption               | 12.544 kWh/24h              |
| Heat emission 100 %              | 828 Watt                    |
| Heat emission default setpoint   | 523 Watt                    |
| Nominal consumption              | 746 Watt                    |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |

## BioPlus EF660W

### General data – BioPlus EF660W

| Technical specifications | Data  |
|--------------------------|---|
| Connection               | 230 VAC, 50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm           |
| Alarm ports              | Voltage-free contact  |
| Access port              | 1 pc. ø24,5 mm  |
| Gross volume             | 660 litres  |
| Net volume               | 484 litres  |
| Door                     | Left or right hinged  |
| Material interior        | Stainless steel and inner doors in PMMA                                   |
| Materiale exterior       | White lacquered steel or stainless steel                                  |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                  |
| Dimensions – W x D x H   | 837 x 756 x 2025/2275 mm  |
| Air system               | BioLine ventilated air distribution system                                |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water              |
| IP class                 | IP21  |



## BioPlus EF660W H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/m²K                  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 360 g                       |
| Refrigeration capacity at -40 °C | 313 Watt                    |
| GWP – CO2e                       | 1411.92                     |
| Energy consumption               | 9.46 kWh/24h                |
| Heat emission 100 %              | 490 Watt                    |
| Heat emission default setpoint   | 394 Watt                    |
| Nominal consumption              | 474 Watt                    |
| Start current                    | 20A                         |
| Sound level                      | –                           |

## BioPlus EF660W H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/m²K                  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 2 x 197 g                   |
| Refrigeration capacity at -40 °C | 2 x 247 Watt                |
| GWP – CO2e                       | 1545.268                    |
| Energy consumption               | 12.92 kWh/24h               |
| Heat emission 100 %              | 963 Watt                    |
| Heat emission default setpoint   | 538 Watt                    |
| Nominal consumption              | 748 Watt                    |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |

## BioPlus EF660W G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/m²K                  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 92 g                        |
| Refrigeration capacity at -40 °C | 244 Watt                    |
| GWP – CO2e                       | N/A                         |
| Energy consumption               | 9.02 kWh/24h                |
| Heat emission 100 %              | 447 Watt                    |
| Heat emission default setpoint   | 376 Watt                    |
| Nominal consumption              | 437 Watt                    |
| Start current                    | 19.5A                       |
| Sound level                      | 55.3 dB(A)                  |

## BioPlus EF660W G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -35/-5 °C                   |
| Ambient temperature range        | +10/+30 °C                  |
| Software variant                 | E1+                         |
| K-Value                          | 0.31 W/m²K                  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 394 g                       |
| Refrigeration capacity at -40 °C | 572 Watt                    |
| GWP – CO2e                       | N/A                         |
| Energy consumption               | 12.87 kWh/24h               |
| Heat emission 100 %              | 839 Watt                    |
| Heat emission default setpoint   | 536 Watt                    |
| Nominal consumption              | 746 Watt                    |
| Start current                    | 14.85A                      |
| Sound level                      | 49 dB(A)                    |

## BioPlus 930

### General data – BioPlus 930

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                        |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24,5 mm   |
| Gross volume             | 930 litres   |
| Net volume               | 702 litres   |
| Door                     | Left or right hinged   |
| Material interior        | Stainless steel  |
| Materiale exterior       | White lacquered steel or stainless steel   |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                               |
| Dimensions – W x D x H   | 780 x 1045 x 2025/2275 mm  |
| Air system               | BioLine ventilated air distribution system   |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                           |
| IP class                 | IP21   |

## BioPlus ER930 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 400 g                       |
| Refrigeration capacity at -10 °C | 487 Watt                    |
| GWP – CO <sub>2</sub> e          | 572                         |
| Energy consumption               | 2.79 kWh/24h                |
| Heat emission 100 %              | 446 Watt                    |
| Heat emission default setpoint   | 116 Watt                    |
| Nominal consumption              | 480 Watt                    |
| Start current                    | 15A                         |
| Sound level                      | –                           |

## BioPlus ER930 H – with glass door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+38 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.509 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 400 g                       |
| Refrigeration capacity at -10 °C | 487 Watt                    |
| GWP – CO <sub>2</sub> e          | 572                         |
| Energy consumption               | 2.80 kWh/24h                |
| Heat emission 100 %              | 423 Watt                    |
| Heat emission default setpoint   | 117 Watt                    |
| Nominal consumption              | 480 Watt                    |
| Start current                    | 15A                         |
| Sound level                      | –                           |



## BioPlus RF930 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 300 g                       |
| Refrigeration capacity at -25 °C | 714 Watt                    |
| GWP – CO <sub>2</sub> e          | 1176.6                      |
| Energy consumption               | 6.28 kWh/24h                |
| Heat emission 100 %              | 668 Watt                    |
| Heat emission default setpoint   | 262 Watt                    |
| Nominal consumption              | 776 Watt                    |
| Start current                    | 18.6A                       |
| Sound level                      | 49.9 dB(A)                  |

## BioPlus ER930 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 126 g                       |
| Refrigeration capacity at -10 °C | 627 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.4158                      |
| Energy consumption               | 2.10 kWh/24h                |
| Heat emission 100 %              | 420 Watt                    |
| Heat emission default setpoint   | 88 Watt                     |
| Nominal consumption              | 457 Watt                    |
| Start current                    | 10.4A                       |
| Sound level                      | –                           |

## BioPlus ER930 G – with glass door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+38° C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.509 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 126 g                       |
| Refrigeration capacity at -10 °C | 627 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.4158                      |
| Energy consumption               | 2.81 kWh/24h                |
| Heat emission 100 %              | 421 Watt                    |
| Heat emission default setpoint   | 117 Watt                    |
| Nominal consumption              | 457 Watt                    |
| Start current                    | 10.4A                       |
| Sound level                      | –                           |





## BioPlus RF930 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T2 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 90 g                        |
| Refrigeration capacity at -25 °C | 828 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.297                       |
| Energy consumption               | 7.25 kWh/24h                |
| Heat emission 100 %              | 675 Watt                    |
| Heat emission default setpoint   | 255 Watt                    |
| Nominal consumption              | 830.5 Watt                  |
| Start current                    | 23.5A                       |
| Sound level                      | –                           |

## BioPlus 1270

### General data – BioPlus 1270

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                        |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24,5 mm   |
| Gross volume             | 1270 litres  |
| Net volume               | 864 litres   |
| Door                     | One left hinged and one right hinged   |
| Material interior        | Stainless steel  |
| Materiale exterior       | White lacquered steel or stainless steel   |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                               |
| Dimensions – W x D x H   | 1390 x 876 x 1875/2125 mm  |
| Air system               | BioLine ventilated air distribution system   |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                           |
| IP class                 | IP21   |



## BioPlus ER1270 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 355 g                       |
| Refrigeration capacity at -10 °C | 703 Watt                    |
| GWP – CO <sub>2</sub> e          | 507.65                      |
| Energy consumption               | 2.85 kWh/24h                |
| Heat emission 100 %              | 578 Watt                    |
| Heat emission default setpoint   | 119 Watt                    |
| Nominal consumption              | 570 Watt                    |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |

## BioPlus ER1270 H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38° C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5752 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 355 g                        |
| Refrigeration capacity at -10 °C | 703 Watt                     |
| GWP – CO <sub>2</sub> e          | 507.65                       |
| Energy consumption               | 3.39 kWh/24h                 |
| Heat emission 100 %              | 578 Watt                     |
| Heat emission default setpoint   | 141 Watt                     |
| Nominal consumption              | 570 Watt                     |
| Start current                    | 14.8A                        |
| Sound level                      | –                            |

---

## BioPlus RF1270 H

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 420 g                       |
| Refrigeration capacity at -25 °C | 947 Watt                    |
| GWP – CO <sub>2</sub> e          | 1647.24                     |
| Energy consumption               | 7.48 kWh/24h                |
| Heat emission 100 %              | 778 Watt                    |
| Heat emission default setpoint   | 312 Watt                    |
| Nominal consumption              | 844 Watt                    |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |



## BioPlus ER1270 H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 2 x 215 g                   |
| Refrigeration capacity at -10 °C | 2 x 703 Watt                |
| GWP – CO <sub>2</sub> e          | 614.9                       |
| Energy consumption               | N/A                         |
| Heat emission 100 %              | N/A                         |
| Heat emission default setpoint   | N/A                         |
| Nominal consumption              | 1040 Watt                   |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |

## BioPlus ER1270 H – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38° C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5752 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 2 x 215 g                    |
| Refrigeration capacity at -10 °C | 2 x 703 Watt                 |
| GWP – CO <sub>2</sub> e          | 614,9                        |
| Energy consumption               | 5.09 kWh/24h                 |
| Heat emission 100 %              | 468 Watt                     |
| Heat emission default setpoint   | 212 Watt                     |
| Nominal consumption              | 1040 Watt                    |
| Start current                    | 14.8A                        |
| Sound level                      | –                            |

## BioPlus RF1270 H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 2 x 220 g                   |
| Refrigeration capacity at -25 °C | 2 x 616 Watt                |
| GWP – CO <sub>2</sub> e          | 1725,68                     |
| Energy consumption               | 8.89 kWh/24h                |
| Heat emission 100 %              | 1048 Watt                   |
| Heat emission default setpoint   | 370 Watt                    |
| Nominal consumption              | 1118 Watt                   |
| Start current                    | 19.6A                       |
| Sound level                      | –                           |



## BioPlus ER1270 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 110 g                       |
| Refrigeration capacity at -10 °C | 946 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.363                       |
| Energy consumption               | 2.50 kWh/24h                |
| Heat emission 100 %              | 530 Watt                    |
| Heat emission default setpoint   | 104 Watt                    |
| Nominal consumption              | 586 Watt                    |
| Start current                    | 13.2A                       |
| Sound level                      | 48.3 dB(A)                  |

## BioPlus ER1270 G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38° C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5752 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 110 g                        |
| Refrigeration capacity at -10 °C | 946 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.363                        |
| Energy consumption               | 3.48 kWh/24h                 |
| Heat emission 100 %              | 540 Watt                     |
| Heat emission default setpoint   | 145 Watt                     |
| Nominal consumption              | 586 Watt                     |
| Start current                    | 13.2A                        |
| Sound level                      | 48.3 dB(A)                   |

## BioPlus RF1270 G

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 85 g                    |
| Refrigeration capacity at -25 °C | 2 x 475 Watt                |
| GWP – CO <sub>2</sub> e          | 0.561                       |
| Energy consumption               | 7.97 kWh/24h                |
| Heat emission 100 %              | 828 Watt                    |
| Heat emission default setpoint   | 332 Watt                    |
| Nominal consumption              | 860 Watt                    |
| Start current                    | 13.2A                       |
| Sound level                      | 51 dB(A)                    |





## BioPlus ER1270 G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 120 g                   |
| Refrigeration capacity at -10 °C | 2 x 719 Watt                |
| GWP – CO <sub>2</sub> e          | 0.792                       |
| Energy consumption               | 2.61 kWh/24h                |
| Heat emission 100 %              | 753 Watt                    |
| Heat emission default setpoint   | 109 Watt                    |
| Nominal consumption              | 828 Watt                    |
| Start current                    | 13.2A                       |
| Sound level                      | 48.3 dB(A)                  |

## BioPlus ER1270 G – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5752 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 2 x 120 g                    |
| Refrigeration capacity at -10 °C | 2 x 719 Watt                 |
| GWP – CO <sub>2</sub> e          | 0.792                        |
| Energy consumption               | 2,945 kWh/24h                |
| Heat emission 100 %              | 751 Watt                     |
| Heat emission default setpoint   | 123 Watt                     |
| Nominal consumption              | 828 Watt                     |
| Start current                    | 13.2A                        |
| Sound level                      | 48.3 dB(A)                   |

## BioPlus RF1270 G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 120 g                   |
| Refrigeration capacity at -25 °C | 2 x 512 Watt                |
| GWP – CO <sub>2</sub> e          | 0.792                       |
| Energy consumption               | 7.974 kWh/24h               |
| Heat emission 100 %              | 828 Watt                    |
| Heat emission default setpoint   | 332 Watt                    |
| Nominal consumption              | 930 Watt                    |
| Start current                    | 13.1A                       |
| Sound level                      | –                           |

## BioPlus 1400

### General data – BioPlus 1400

| Technical specifications | Data   |
|--------------------------|--|
| Connection               | 230 VAC/50 Hz  |
| Control Unit             | Gram Control Unit with voltage-free contact, E-sensor, dry cooling and offset function |
| Alarms                   | Acoustic and visual, High/Low temperature alarms and door alarm                        |
| Alarm ports              | Voltage-free contact   |
| Access port              | 1 pc. ø24,5 mm   |
| Gross volume             | 1400 litres  |
| Net volume               | 968 litres   |
| Door                     | One left hinged and one right hinged   |
| Material interior        | Stainless steel  |
| Materiale exterior       | White lacquered steel or stainless steel   |
| Insulation               | 60 mm polyurethane with HFC-free cyclopentane propellant                               |
| Dimensions – W x D x H   | 1390 x 876 x 2025/2275 mm  |
| Air system               | BioLine ventilated air distribution system   |
| Defrost system           | Automatic smart defrost with re-evaporation of defrost water                           |
| IP class                 | IP21   |

## BioPlus ER1400 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 355 g                       |
| Refrigeration capacity at -10 °C | 703 Watt                    |
| GWP – CO <sub>2</sub> e          | 507,65                      |
| Energy consumption               | 2.80 kWh/24h                |
| Heat emission 100 %              | 553 Watt                    |
| Heat emission default setpoint   | 117 Watt                    |
| Nominal consumption              | 570 Watt                    |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |

## BioPlus ER1400 H – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38° C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5808 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R134a                        |
| Refrigerant charge               | 355 g                        |
| Refrigeration capacity at -10 °C | 703 Watt                     |
| GWP – CO <sub>2</sub> e          | 507,65                       |
| Energy consumption               | 4.16kWh/24h                  |
| Heat emission 100 %              | 577 Watt                     |
| Heat emission default setpoint   | 173 Watt                     |
| Nominal consumption              | 570 Watt                     |
| Start current                    | 14.8A                        |
| Sound level                      | –                            |



## BioPlus RF1400 H – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 420 g                       |
| Refrigeration capacity at -25 °C | 947 Watt                    |
| GWP – CO <sub>2</sub> e          | 1647,24                     |
| Energy consumption               | 8.08 kWh/24h                |
| Heat emission 100 %              | 815 Watt                    |
| Heat emission default setpoint   | 337 Watt                    |
| Nominal consumption              | 844 Watt                    |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |



## BioPlus ER1400 H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 2 x 215 g                   |
| Refrigeration capacity at -10 °C | 2 x 703 Watt                |
| GWP – CO <sub>2</sub> e          | 614.9                       |
| Energy consumption               | 2.83 kWh/24h                |
| Heat emission 100 %              | 462 Watt                    |
| Heat emission default setpoint   | 118 Watt                    |
| Nominal consumption              | 1040 Watt                   |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |

## BioPlus ER1400 H – with dual compressor and glass door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+38° C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.5808 W/(m2*K)             |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R134a                       |
| Refrigerant charge               | 2 x 215 g                   |
| Refrigeration capacity at -10 °C | 2 x 703 Watt                |
| GWP – CO2e                       | 614.9                       |
| Energy consumption               | 3.90 kWh/24h                |
| Heat emission 100 %              | 452 Watt                    |
| Heat emission default setpoint   | 162 Watt                    |
| Nominal consumption              | 1040 Watt                   |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |

## BioPlus RF1400 H – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m2*K)               |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R404A                       |
| Refrigerant charge               | 2 x 220 g                   |
| Refrigeration capacity at -25 °C | 2 x 616 Watt                |
| GWP – CO2e                       | 1725.68                     |
| Energy consumption               | 9.45 kWh/24h                |
| Heat emission 100 %              | 1071 Watt                   |
| Heat emission default setpoint   | 394 Watt                    |
| Nominal consumption              | 1118 Watt                   |
| Start current                    | 14.8A                       |
| Sound level                      | –                           |



## BioPlus ER1400 G – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 110 g                       |
| Refrigeration capacity at -10 °C | 946 Watt                    |
| GWP – CO <sub>2</sub> e          | 0.363                       |
| Energy consumption               | 2.56 kWh/24h                |
| Heat emission 100 %              | 538 Watt                    |
| Heat emission default setpoint   | 106 Watt                    |
| Nominal consumption              | 586 Watt                    |
| Start current                    | 13.2A                       |
| Sound level                      | 47.5 dB(A)                  |

## BioPlus ER1400 G – with glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38° C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5808 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 110 g                        |
| Refrigeration capacity at -10 °C | 946 Watt                     |
| GWP – CO <sub>2</sub> e          | 0.363                        |
| Energy consumption               | 3.72 kWh/24h                 |
| Heat emission 100 %              | 539 Watt                     |
| Heat emission default setpoint   | 155 Watt                     |
| Nominal consumption              | 586 Watt                     |
| Start current                    | 13.2A                        |
| Sound level                      | 47.5 dB(A)                   |

## BioPlus RF1400 G

### – with solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 85 g                    |
| Refrigeration capacity at -25 °C | 2 x 475 Watt                |
| GWP – CO <sub>2</sub> e          | 0.561                       |
| Energy consumption               | 8.39 kWh/24h                |
| Heat emission 100 %              | 790 Watt                    |
| Heat emission default setpoint   | 349 Watt                    |
| Nominal consumption              | 860 Watt                    |
| Start current                    | 13.2A                       |
| Sound level                      | 46.6 dB(A)                  |





## BioPlus ER1400 G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -2/+20 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | M5+                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 120 g                   |
| Refrigeration capacity at -10 °C | 2 x 719 Watt                |
| GWP – CO <sub>2</sub> e          | 0.792                       |
| Energy consumption               | 2.81 kWh/24h                |
| Heat emission 100 %              | 761 Watt                    |
| Heat emission default setpoint   | 117 Watt                    |
| Nominal consumption              | 828 Watt                    |
| Start current                    | 13.6A                       |
| Sound level                      | –                           |



## BioPlus ER1400 G – with dual compressor and glass door

| Technical specifications         | Data                         |
|----------------------------------|------------------------------|
| Temperature range                | -2/+20 °C                    |
| Ambient temperature range        | +10/+38 °C                   |
| Software variant                 | M5+                          |
| K-Value                          | 0.5808 W/(m <sup>2</sup> *K) |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc  |
| ATEX certificate                 | –                            |
| Refrigerant                      | R290                         |
| Refrigerant charge               | 2 x 120                      |
| Refrigeration capacity at -10 °C | 2 x 719 Watt                 |
| GWP – CO <sub>2</sub> e          | 0.792                        |
| Energy consumption               | 3.965 kWh/24h                |
| Heat emission 100 %              | 767 Watt                     |
| Heat emission default setpoint   | 165 Watt                     |
| Nominal consumption              | 828 Watt                     |
| Start current                    | 13.6A                        |
| Sound level                      | –                            |

## BioPlus RF1400 G – with dual compressor and solid door

| Technical specifications         | Data                        |
|----------------------------------|-----------------------------|
| Temperature range                | -25/-5 °C                   |
| Ambient temperature range        | +10/+43 °C                  |
| Software variant                 | F51                         |
| K-Value                          | 0.31 W/(m <sup>2</sup> *K)  |
| ATEX Marking                     | II 3G Ex nA nC nL IIB T3 Gc |
| ATEX certificate                 | –                           |
| Refrigerant                      | R290                        |
| Refrigerant charge               | 2 x 85 g                    |
| Refrigeration capacity at -25 °C | 2 x 633 Watt                |
| GWP – CO <sub>2</sub> e          | 0.561                       |
| Energy consumption               | 8.385 kWh/24h               |
| Heat emission 100 %              | 790 Watt                    |
| Heat emission default setpoint   | 349 Watt                    |
| Nominal consumption              | 930 Watt                    |
| Start current                    | 13.2A                       |
| Sound level                      | 46.6 dB(A)                  |



# Declaration of conformity

## BioMidi



### English EC Declaration of Conformity

We, **Gram Scientific ApS**, declare as manufacturers under sole responsibility that the following products comply with all relevant regulations:

|                         |   |
|-------------------------|---|
| <b>Range:</b>           | <b>BioMidi</b>  |
| Model:                  | RR425, RF425, RR625, RF625 & EF425                      |
| Refrigeration:          | R290, R404A & R134a                                     |
| Product description:    | Refrigerators and freezer for laboratory and biostorage |
| Valid from (Year/Week): | 2023/01   |

This declaration pertains to compliance with all applicable essential requirements and other provisions of the European Council Directive and regulations. Specifically, the following Directives and Regulations of the European Parliament and of the Council apply:

| Directives and Regulations of the European Parliament and of the Council: |  |
|---|--|
| Machinery Directive 2006/42/EC  |  |
| - ATEX Directive 2014/34/EU   |  |
| - Pressure Equipment Directive 2014/68/EU                                 |  |
| - Low Voltage Directive 2014/35/EU  |  |
| - EMC Directive 2014/30/EU  |  |
| - RoHS Directive 2011/65/EU   |  |
| - REACH EC No.1907/2006   |  |
| - F-Gas Regulation (EU) No 2024/573                                       |  |

Product compliance has been demonstrated based on the following harmonized standards:

| Harmonized Standards: | Text:  |
|-----------------------|--|
| EN 61010-1:2010       | Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements   |
| EN 60601-1:2006       | Medical electrical equipment. General requirements for basic safety and essential performance  |
| EN 60601-1-2:2015     | Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests |
| EN 60079-0:2012       | Electrical apparatus for explosive atmospheres – Part 0: General requirements  |
| EN 60079-11:2012      | Explosive atmospheres – Part 11: Equipment protection by intrinsic safety "I"  |
| EN 60079-15:2010      | Explosive atmospheres – Part 15: Equipment protection by type of protection "n"  |
| EN 60079-25:2010      | Explosive atmospheres – Part 25: Intrinsically safe systems  |
| EN ISO 3744:2010      | Acoustics – Determination of sound power levels of noise sources using sound pressure - Engineering method in an essentially free field over a reflecting plane                        |
| EN ISO 9001:2015      | Quality management systems   |
| EN ISO 14001:2015     | Environment management systems – Requirements with guidance for use  |

Gram Scientific ApS  
Aage Grams Vej 1  
DK-6500 Vejens  
Denmark  
Telephone: + 45 73 20 13 00

Vejens, 05.03.2024

John B. S. Petersen  
Approval Manager

# BioMidi – Accessory code 69



## English EC Declaration of Conformity

We, **Gram Scientific ApS**, declare as manufacturer under sole responsibility that the following products comply with all relevant regulations:

|                         |   |
|-------------------------|---|
| <b>Range:</b>           | <b>BioMidi (Accessory code 69)</b>          |
| Model:                  | RR425, RF425, RR625, RF625 & EF425          |
| Refrigeration:          | R290, R404A & R134a                         |
| Product description:    | Refrigerators for laboratory and biostorage |
| Valid from (Year/Week): | 2024/01                                     |

This declaration pertains to compliance with all applicable essential requirements and other provisions of the European Council Directive and regulations. Specifically, the following Directives and Regulations of the European Parliament and of the Council apply:

**Directives and Regulations of the European Parliament and of the Council:**

Machinery Directive 2006/42/EC  
 - Pressure Equipment Directive 2014/68/EU  
 - Low Voltage Directive 2014/35/EU  
 - EMC Directive 2014/30/EU  
 - RoHS Directive 2011/65/EU  
 - REACH EC No.1907/2006  
 - F-Gas Regulation (EU) No 2024/573

Product compliance has been demonstrated based on the following harmonized standards:

| Harmonized Standards: | Text:  |
|-----------------------|--|
| EN 61010-1:2010       | Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements   |
| EN 60601-1:2006       | Medical electrical equipment. General requirements for basic safety and essential performance  |
| EN 60601-1-2:2015     | Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests |
| DIN 13277:2022-05     | Refrigerators and freezers for laboratory and medical applications – Terminology, requirements, testing  |
| EN ISO 3744:2010      | Acoustics – Determination of sound power levels of noise sources using sound pressure - Engineering method in an essentially free field over a reflecting plane                        |
| EN ISO 9001:2015      | Quality management systems   |
| EN ISO 14001:2015     | Environment management systems – Requirements with guidance for use  |

Gram Scientific ApS  
 Aage Grams Vej 1  
 DK-6500 Vojens  
 Denmark  
 Telephone: + 45 73 20 13 00

Vojens, 15.05.2024



John B. S. Petersen  
 Approval Manager



## English EC Declaration of Conformity

We, **Gram Scientific ApS**, declare as manufacturers under sole responsibility that the following products comply with all relevant regulations:

|                                |   |
|--------------------------------|---|
| <b>Range:</b>                  | <b>BioPlus</b>  |
| <b>Model:</b>                  | ER500, RF500, ER600D, RF600D, ER600W, RF600W, ER660D, RF660D, ER660W, RF660W, ER930, RF930, ER1270, RF1270, ER1400, RF1400, EF600W & EF660W |
| <b>Refrigeration:</b>          | R290, R404A & R134a   |
| <b>Product description:</b>    | Refrigerators and freezer for laboratory and biostorage   |
| <b>Valid from (Year/Week):</b> | 2023/01   |

This declaration pertains to compliance with all applicable essential requirements and other provisions of the European Council Directive and regulations. Specifically, the following Directives and Regulations of the European Parliament and of the Council apply:

### Directives and Regulations of the European Parliament and of the Council:

Machinery Directive 2006/42/EC  
 - ATEX Directive 2014/34/EU  
 - Pressure Equipment Directive 2014/68/EU  
 - Low Voltage Directive 2014/35/EU  
 - EMC Directive 2014/30/EU  
 - RoHS Directive 2011/65/EU  
 - REACH EC No.1907/2006  
 - F-Gas Regulation (EU) No 2024/573

Product compliance has been demonstrated based on the following harmonized standards:

| Harmonized Standards: | Text:  |
|-----------------------|--|
| EN 61010-1:2010       | Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements   |
| EN 60601-1:2006       | Medical electrical equipment. General requirements for basic safety and essential performance  |
| EN 60601-1-2:2015     | Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests |
| EN 60079-0:2012       | Electrical apparatus for explosive atmospheres – Part 0: General requirements  |
| EN 60079-11:2012      | Explosive atmospheres – Part 11: Equipment protection by intrinsic safety "I"  |
| EN 60079-15:2010      | Explosive atmospheres – Part 15: Equipment protection by type of protection "n"  |
| EN 60079-25:2010      | Explosive atmospheres – Part 25: Intrinsically safe systems  |
| EN ISO 3744:2010      | Acoustics – Determination of sound power levels of noise sources using sound pressure - Engineering method in an essentially free field over a reflecting plane                        |
| EN ISO 9001:2015      | Quality management systems   |
| EN ISO 14001:2015     | Environment management systems – Requirements with guidance for use  |

Gram Scientific ApS  
 Aage Grams Vej 1  
 DK-6500 Vejens  
 Denmark  
 Telephone: + 45 73 20 13 00

Vejens, 05.03.2024

John B. S. Petersen  
 Approval Manager

# BioPlus – Accessory code 69



## English EC Declaration of Conformity

We, **Gram Scientific ApS**, declare as manufacturer under sole responsibility that the following products comply with all relevant regulations:

|                                |   |
|--------------------------------|---|
| <b>Range:</b>                  | <b>BioPlus (Accessory code 69)</b>  |
| <b>Model:</b>                  | ER500, RF500, ER600D, RF600D, ER600W, RF600W, ER660D, RF660D, ER660W, RF660W, ER930, RF930, ER1270, RF1270, ER1400, RF1400, EF600W & EF660W |
| <b>Refrigeration:</b>          | R290, R404A & R134a   |
| <b>Product description:</b>    | Refrigerators for laboratory and biostorage   |
| <b>Valid from (Year/Week):</b> | 2024/01   |

This declaration pertains to compliance with all applicable essential requirements and other provisions of the European Council Directive and regulations. Specifically, the following Directives and Regulations of the European Parliament and of the Council apply:

| Directives and Regulations of the European Parliament and of the Council: |  |
|---|--|
| Machinery Directive 2006/42/EC  |  |
| - Pressure Equipment Directive 2014/68/EU                                 |  |
| - Low Voltage Directive 2014/35/EU  |  |
| - EMC Directive 2014/30/EU  |  |
| - RoHS Directive 2011/65/EU   |  |
| - REACH EC No.1907/2006   |  |
| - F-Gas Regulation (EU) No 2024/573                                       |  |

Product compliance has been demonstrated based on the following harmonized standards:

| Harmonized Standards: | Text:  |
|-----------------------|--|
| EN 61010-1:2010       | Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements   |
| EN 60601-1:2006       | Medical electrical equipment. General requirements for basic safety and essential performance  |
| EN 60601-1-2:2015     | Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests |
| DIN 13277:2022-05     | Refrigerators and freezers for laboratory and medical applications – Terminology, requirements, testing  |
| EN ISO 3744:2010      | Acoustics – Determination of sound power levels of noise sources using sound pressure - Engineering method in an essentially free field over a reflecting plane                        |
| EN ISO 9001:2015      | Quality management systems   |
| EN ISO 14001:2015     | Environment management systems – Requirements with guidance for use  |


Gram Scientific ApS  
 Aage Grams Vej 1  
 DK-6500 Vejens  
 Denmark  
 Telephone: + 45 73 20 13 00

Vejens, 15.05.2024


  
 John B. S. Petersen  
 Approval Manager

# Piping diagram

## BioMidi/BioPlus




|   | DK             | GB             | D                |
|---|----------------|----------------|------------------|
| A | Kompressor     | Compressor     | Kompressor       |
| B | Kondensator    | Condenser      | Verflüssiger     |
| C | Tørrefilter    | Filter drier   | Trockenfilter    |
| D | Varmeudveksler | Heat exchanger | Wärmeaustauscher |
| E | Fordamper      | Evaporator     | Verdampfer       |

|   |              |   |            |      |          |
|---|--------------|---|------------|------|----------|
| Description: Rørdiagram/Piping diagram/Rohrleitungsplan HAV   |              | General tolerances:<br>ISO 2768-1 (m)<br>ISO 2768-2 (L) |            |      |          |
| <br>© Gram Scientific ApS<br>DENMARK | Date: 141201 | Drawing no.:  | Sheet no.: | Rev. | A format |
|   | Name: JP     | 765040844   | 1<br>of 1  | 00   | 4        |
|   | Scale 1/1    |   |            |      |          |

All rights to this drawing belongs to Gram Scientific ApS, cf. law of Copyright. Thus the drawing should not without our written permission be copied, presented or passed over to a third person. Misuse will lead to prosecu




## BioPlus – With dual refrigeration




|   | DK             | GB             | D                |
|---|----------------|----------------|------------------|
| A | Kompressor     | Compressor     | Kompressor       |
| B | Kondensator    | Condenser      | Verflüssiger     |
| C | Tørrefilter    | Filter drier   | Trockenfilter    |
| D | Varmeudveksler | Heat exchanger | Wärmeaustauscher |
| E | Fordamper      | Evaporator     | Verdampfer       |

Revision:


|   |   |              |   |      |        |
|---|---|--------------|---|------|--------|
| Description: Piping diagram BioPlus w/dual refrigeration  |   |              | General tolerances:<br>ISO 2768-1 (m)<br>ISO 2768-2 (L) |      |        |
| <br>© Gram Scientific ApS<br>DENMARK | Date: 09-07-2010  | Drawing no.: | Sheet no.:  | Rev. | A form |
|   | Name: JABP  | 765041375    | 1<br>of<br>1  | 000  | 4      |
|   | Scale 1/1   |              |   |      |        |
|   | All rights to this drawing belongs to Gram Scientific ApS, cf. law of Copyright. Thus the drawing should not without our written permission be copied, presented or passed over to a third person. Misuse will lead to prosecution. |              |   |      |        |

# Wiring diagram


## BioPlus RF – With SSR



## BioPlus ER – With SSR



# BioPlus EF600/660 – With SSR



Revision:

Description: Ledn.digram EF/PF 600/660 m/SSR


General Information:  
ISO 2768-1 (M)  
ISO 2768-2 (L)

| Doc. no.   | Rev. | Sheet no. | Rev. |
|------------|------|-----------|------|
| 25-01-2016 | 1/1  | 765072826 | 1/1  |


© Gram Scientific ApS  
DENMARK

All rights to this drawing belong to Gram Scientific ApS. All use of copyright must be the drawing holder's responsibility. No part of this drawing may be reproduced or transmitted in any form or by any means electronic, mechanical, photocopying, recording, or by any information storage and retrieval system, without permission in writing from Gram Scientific ApS.

## BioPlus ER1270/1400 – With SSR




# BioPlus RF1270/1400 - With dual compressor - With SSR




| Revision   |  |               |  |
|--|--|---------------|--|
| Description: Lehn, diagram RF 1270/1400 2 kompr. m/SSR |  |               |  |
| Drawing no.:   |  | Drawing no.:  |  |
| Date: 10-2016  |  | Date: 10-2016 |  |
| Scale: 1/1   |  | Scale: 1/1    |  |
| Sheet no.:   |  | Sheet no.:    |  |
| 1  |  | 1             |  |
| 1000   |  | 1000          |  |
| 2  |  | 2             |  |

All rights to this drawing belongs to Gram Scientific A/S. All law of copyright. This drawing should not be used without the written permission of Gram Scientific A/S. The drawing is for information only. It is not a contract. The drawing is for information only. It is not a contract. The drawing is for information only. It is not a contract.

# BioPlus ER1270/1400 - With dual compressor - With SSR




# BioPlus RF/EF - With dual compressor - With SSR






# BioPlus ER1270/1400 - With dual compressor - With LTP and SSR




# BioMidi RR425/625 – With LTP




# BioMidi RF425/625 – With SSR

Wire colors:  
L = brown  
N = blue  
W = white  
R = red  
B = black



|  |  |  |  |
|--|--|--|--|
| Revision:  |  |  |  |
| Description: Ledn:diagram RF 425/625 G/H m/SSR   |  |  |  |
| General dimensions   |  |  |  |
| ISO 2788-1 (n)   |  |  |  |
| ISO 2788-2 (l)   |  |  |  |
| Date: 02-05-2016   |  |  |  |
| Scale: 1/1   |  |  |  |
| Drawing no.: 765072833   |  |  |  |
| Sheet no.: 1   |  |  |  |
| Total: 1   |  |  |  |
| Author: 2  |  |  |  |
| Check: 1   |  |  |  |
| Drawn: 1   |  |  |  |
| Approved: 1  |  |  |  |
| All rights to this drawing belong to Gram Scientific Aps. If law of copyright, thus the drawing should not be reproduced without written permission by Gram Scientific Aps. If the drawing is used for reproduction, the user will need to pay for reproduction. |  |  |  |

## 152





## Installation Qualification Operation Qualification

The following IQ/OQ is intended to be a guideline, local IQ/OQ procedures can vary depending on application and items stored in the Gram BioLine cabinet.

Deviations from the specifications dictated in the PQ are to be reported in the deviation report.

The IQ/OQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

Organisation:

Location of installation:

Model:

Serial number:

Item and revision number of instructions for use:

Status of operation:

- ☐ Active
- ☐ Inactive

Name of vendor:

Warranty:

Start:

End:

Instructions on use to starting the cabinet:

1. Training of the responsible party Date: \_\_\_\_\_ By: \_\_\_\_\_

2. Operational test of the cabinet Date: \_\_\_\_\_ By: \_\_\_\_\_

3. Responsible party \_\_\_\_\_ Tel: \_\_\_\_\_

Instructions to users:

*The responsible party is trained in use of the cabinet in reference to the user manual*

☐ General use of cabinet

Objections to the mentioned:

☐ Service & maintenance

☐ The cabinet was delivered without defects/damage.  
The cabinet started as specified in the user manual

---

---

---

---

#### Set values:

☐ Setpoint temperature \_\_\_\_\_ °C

#### Local alarm settings:

☐ Upper temperature alarm limit (LhL) \_\_\_\_\_ °C

☐ Lower temperature alarm limit (LLL) \_\_\_\_\_ °C

#### External alarm settings:

(See voltage-free contact in user manual)

☐ Upper temperature alarm limit (EhL) \_\_\_\_\_ °C

☐ Lower temperature alarm limit (ELL) \_\_\_\_\_ °C

#### Factory settings

| Model/<br>Setpoint temperature |        | LhL    | LLL    | EhL    | ELL    |
|--------------------------------|--------|--------|--------|--------|--------|
| RR – BioMidi                   | +5 °C  | +25 °C | 0 °C   | +25 °C | 0 °C   |
| ER – BioPlus                   | +5 °C  | +25 °C | -5 °C  | +25 °C | -5 °C  |
| RF – BioPlus/BioMidi           | -20 °C | +25 °C | -35 °C | +25 °C | -35 °C |
| EF – BioPlus                   | -35 °C | +25 °C | -45 °C | +25 °C | -45 °C |
| EF – BioMidi                   | -40 °C | +25 °C | -60 °C | +25 °C | -60 °C |

|       |                       |            |                     |            |
|-------|-----------------------|------------|---------------------|------------|
| Date: | Name of trained user: | Signature: | Name of instructor: | Signature: |
| _____ | _____                 | _____      | _____               | _____      |
| _____ | _____                 | _____      | _____               | _____      |
| _____ | _____                 | _____      | _____               | _____      |
| _____ | _____                 | _____      | _____               | _____      |

Model: \_\_\_\_\_ SN: \_\_\_\_\_

## Installation Qualification – IQ

| ID   | Description of installation   | Reference in manual | Comply |    | Attachment | Notes |
|------|---|---------------------|--------|----|------------|-------|
|      |   |                     | Yes    | No |            |       |
| I-1  | Ensure the cabinet is installed in-doors.   | page 10             |        |    |            |       |
| I-2  | Ensure the cabinet is installed in a dry and sufficiently ventilated area.                  | page 10             |        |    |            |       |
| I-3  | Ensure the cabinet is not in direct contact with sunlight or other heat sources.            | page 10             |        |    |            |       |
| I-4  | Ensure that the ambient operating temperature is within the allowed range.                  | page 10             |        |    |            |       |
| I-5  | Ensure that the cabinet is not installed in a chloric/acidic environment.                   | page 10             |        |    |            |       |
| I-6  | Ensure that the protective film on the cabinet is removed.                                  | page 10             |        |    |            |       |
| I-7  | Ensure that the cabinet is cleaned with a mild soap solution                                | page 10             |        |    |            |       |
| I-8  | Ensure that the cabinet has stood upright for 24 hours if the cabinet has been laying down. | page 10             |        |    |            |       |
| I-9  | Ensure that the cabinet is levelled if it is equipped with legs.                            | page 11             |        |    |            |       |
| I-10 | Ensure a level surface if the cabinet is equipped with wheels/castors.                      | page 11             |        |    |            |       |
| I-11 | If equipped with wheels/castors: Ensure wheels/castors are locked after positioning         | page 11             |        |    |            |       |
| I-12 | If equipped with drawers and/or glass door: Ensure that tilt-bracket is mounted.            | page 12             |        |    |            |       |
| I-13 | Ensure a distance of 15-75 mm between cabinet and back wall.                                | page 13             |        |    |            |       |
| I-14 | Ensure that there is minimum a gap of 30 mm between cabinets.                               | page 14             |        |    |            |       |
| I-15 | Ensure that the upper part of the cabinet is not covered.                                   | page 14             |        |    |            |       |
| I-16 | Ensure that electrical appliances are not being used inside the cabinet.                    | page 14             |        |    |            |       |
| I-17 | Ensure connection from voltage-free contact to external monitoring system                   | page 15             |        |    |            |       |

Model: \_\_\_\_\_

SN: \_\_\_\_\_



## Installation Qualification - IQ

| ID     | Description of installation  | Reference in manual | Comply |    | Attachment | Notes |
|--------|--|---------------------|--------|----|------------|-------|
|        |  |                     | Yes    | No |            |       |
| I-18   | Ensure that the inner doors can operate in accordance with the instructions.           | N/A                 |        |    |            |       |
| I-19   | Ensure the correct electrical connection (compare local values with type/number plate) | page 16             |        |    |            |       |
| I-20-1 | Ensure that the power cord is secured by the preload cover                             | page 16             |        |    |            |       |
| I-20-2 | Ensure that the power cord is secured by the preload hanger                            | N/A                 |        |    |            |       |
| I-21   | Mark the power cord with: "Do not separate when energized".                            | page 16             |        |    |            |       |
| I-22-1 | Ensure equipotential bonding (applicable for ATEX Cat. 3 Zone 2 areas)                 | page 18             |        |    |            |       |

Model: \_\_\_\_\_

SN: \_\_\_\_\_

## Operation Qualification – OQ

| ID   | Description of operation   | Reference in manual | Comply |    | Attachment | Notes |
|------|--|---------------------|--------|----|------------|-------|
|      |  |                     | Yes    | No |            |       |
| O-1  | Turn on the cabinet – Display test (software version and variant). | page 34             |        |    |            |       |
| O-2  | Set/adjust set-point temperature.                                  | page 34             |        |    |            |       |
| O-3  | Set/adjust LhL – Upper alarm limit (local).                        | page 38             |        |    |            |       |
| O-4  | Set/adjust LLL – Lower alarm limit (local).                        | page 38             |        |    |            |       |
| O-5  | Set/adjust Lhd – Delay of the upper alarm limit (local).           | page 39             |        |    |            |       |
| O-6  | Set/adjust LLd – Delay of the lower alarm limit (local).           | page 39             |        |    |            |       |
| O-7  | Activate / deactivate dA – Door alarm (local).                     | page 40             |        |    |            |       |
| O-8  | Set/adjust dAd – Delay of the door alarm (local).                  | page 40             |        |    |            |       |
| O-9  | Activate / deactivate BU – Acoustic alarms (local).                | page 41             |        |    |            |       |
| O-10 | Set/adjust EhL – Upper alarm limit (external).                     | page 42             |        |    |            |       |
| O-11 | Set/adjust ELL – Lower alarm limit (external).                     | page 42             |        |    |            |       |
| O-12 | Set/adjust Ehd – Delay of the upper alarm limit (external).        | page 43             |        |    |            |       |
| O-13 | Set/adjust ELd – Delay of the lower alarm limit (external).        | page 43             |        |    |            |       |
| O-14 | Activate / deactivate dA – Door alarm (external)                   | page 44             |        |    |            |       |
| O-15 | Set/adjust dAd – Delay of the door alarm (external).               | page 44             |        |    |            |       |
| O-16 | Activate / deactivate BU – Acoustic external alarms.               | page 45             |        |    |            |       |
| O-17 | Set/adjust defrost cycles (dEF) per 24 hours (factory setting: 4). | page 49             |        |    |            |       |
| O-18 | Select reference sensor for the display (dPS) (A or E).            | page 50             |        |    |            |       |

Model: \_\_\_\_\_

SN: \_\_\_\_\_



---

Approval of the test results – Installation Qualification (IQ)

- ☐ The steps in the Installation Qualification – IQ were completed with positive results
- ☐ The steps in the Installation Qualification – IQ were completed with negative results

ID of steps with negative results: \_\_\_\_\_

Approval of the test results – Operation Qualification (OQ)

- ☐ The steps in the Operation Qualification – OQ were completed with positive results
- ☐ The steps in the Operation Qualification – OQ were completed with negative results

ID of steps with negative results: \_\_\_\_\_

Organisation/Responsible party: \_\_\_\_\_

Trainer/Responsible party: \_\_\_\_\_

\_\_\_\_\_  
Stamp & Signature

\_\_\_\_\_  
Stamp & Signature

\_\_\_\_\_  
Tel.

\_\_\_\_\_  
Tel.

\_\_\_\_\_  
e-mail

\_\_\_\_\_  
e-mail

\_\_\_\_\_  
Location & Date

\_\_\_\_\_  
Location & Date

Model: \_\_\_\_\_ SN: \_\_\_\_\_



NOTES:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

## Performance Qualification

Organisation:

---

---

---

Location of installation:

---

---

---

Model: \_\_\_\_\_

SN: \_\_\_\_\_

Item number: \_\_\_\_\_  
(manual)

The PQ consists of inspections of the correct operation of the cabinet under predefined conditions and procedures.

Prerequisites for the PQ are IQ (Installation Qualification) and OQ (Operation Qualification), these must be concluded successfully prior to the initiation of the PQ.

Person responsible for the cabinet:

Name: \_\_\_\_\_  
Date: \_\_\_\_\_  
Signature: \_\_\_\_\_

Person responsible for the test:

Name: \_\_\_\_\_  
Date: \_\_\_\_\_  
Organisation: \_\_\_\_\_  
Signature: \_\_\_\_\_

Person responsible for verification of the test:

Name: \_\_\_\_\_  
Date: \_\_\_\_\_  
Organisation: \_\_\_\_\_  
Signature: \_\_\_\_\_

Test duration:

Initiation (date/time): \_\_\_\_\_  
Conclusion (date/time): \_\_\_\_\_  
Date: \_\_\_\_\_  
Organisation: \_\_\_\_\_  
Signature: \_\_\_\_\_

Model: \_\_\_\_\_

SN: \_\_\_\_\_

List of names – Persons involved in the test procedure and subsequent report

[illegible]

Model: \_\_\_\_\_ SN: \_\_\_\_\_

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.  
The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

| Measurements – Prerequisites |  |     |    |
|------------------------------|--|-----|----|
|                              |  | Yes | No |
| P-1                          | <p>The cabinet must be empty while conducting tests, ie without interior fittings such as drawers, shelves etc.</p> <p>Attachment:</p> <p>Notes:</p> |     |    |
| P-2                          | <p>The measurements must be conducted in accordance to IEC 60068-3-5.</p> <p>Attachment:</p> <p>Notes:</p>   |     |    |
| P-3                          | <p>The positioning of the sensors in the cabinet must be documented with a sketch and/or a photograph.</p> <p>Attachment:</p> <p>Notes:</p>          |     |    |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_  
 Conducted by: \_\_\_\_\_  
 Inspected/verified by: \_\_\_\_\_  
 Model: \_\_\_\_\_ SN: \_\_\_\_\_



Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.  
The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

| Measurements – Prerequisites |   |     |    |
|------------------------------|---|-----|----|
|                              |   | Yes | No |
| P-4                          | <p>Measurements made during the PQ tests must be documented and attached to the PQ.</p> <p>Attachment:</p> <p>Notes:</p>  |     |    |
| P-5                          | <p>Specify setpoint temperature: _____ °C</p> <p>Specify the ambient temperature: _____ °C</p> <p>Attachment:</p> <p>Notes:</p>   |     |    |
| P-6                          | <p>Allowed tolerances –<br/><i>Select the tolerance, according to the model being tested.</i><br/><i>Find model-specific tolerances in appendix.</i></p> <p>Tolerance: +/- _____ K</p> <p>Attachment:</p> <p>Notes:</p> |     |    |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_  
 Conducted by: \_\_\_\_\_  
 Inspected/verified by: \_\_\_\_\_  
 Model: \_\_\_\_\_ SN: \_\_\_\_\_

Deviations from the specifications dictated in the PQ are to be reported in the deviation report.  
The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

| Measurements – Temperature stabilisation |  |     |    |
|--|--|-----|----|
|  |  | Yes | No |
| P-7                                      | <p>The test is intended to provide substantiation for the temperature stability inside the cabinet during normal operation.</p> <p>The temperature inside the cabinet must be stabilised – Where all the points in the working space have reached and maintained the same temperature.</p> <p>When the system is stable, document ordinary operation of the cabinet at the setpoint temperature and ambient temperature specified in P-5.</p> <p>Duration: _____</p> <p>The measurements throughout the operation test, must be documented and attached the PQ.</p> <p>Attachment:</p> <p>Notes:</p> |     |    |
| P-8                                      | <p>Are the measurements inside the allowed tolerances specified in P-6 ?</p> <p>Attachment:</p> <p>Notes:</p>  |     |    |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_  
 Conducted by: \_\_\_\_\_  
 Inspected/verified by: \_\_\_\_\_  
 Model: \_\_\_\_\_ SN: \_\_\_\_\_

Deviations from the specifications dictated in the PQ are to be reported in the deviation report.  
The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

| Measurements – Door opening test |  |     |    |
|----------------------------------|--|-----|----|
|                                  |  | Yes | No |
| P-9                              | <p>The test is intended to provide substantiation for the temperature recovery time inside the cabinet subsequently after a door opening.</p> <p>The temperature inside the cabinet must be stabilised – Where all the points in the working space have reached and maintained the same temperature, the setpoint temperature is specified in P-5.</p> <p>When the system is stable, open the door at 90° for 60 seconds.</p> <p>The measurements, throughout the door opening test, must be documented and attached the PQ.</p> <p>Duration: _____</p> <p>Attachment:</p> <p>Notes:</p> |     |    |
| P-10                             | <p>Has the setpoint temperature specified in P-5, measured in the absolute centre of the cabinet, been achieved within the set time-frame specified in the appendix?</p> <p>Attachment:</p> <p>Notes:</p>  |     |    |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_  
 Conducted by: \_\_\_\_\_  
 Inspected/verified by: \_\_\_\_\_  
 Model: \_\_\_\_\_ SN: \_\_\_\_\_

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.  
The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

| Measurements – Pull-down |   |     |    |
|--------------------------|---|-----|----|
|                          |   | Yes | No |
| P-11                     | <p>The test is intended to provide substantiation for the time it takes for the inside of the cabinet to reach the setpoint temperature specified in P-5.<br/>The initial temperature in the working space is the ambient temperature specified in P-5.<br/>The temperature inside the cabinet must be stabilised in all points of the working space.</p> <p>When the system is stable, turn on the power to the cabinet.</p> <p>The measurements, throughout the pull-down test, must be documented and attached the PQ.</p> <p>Duration: _____</p> <p>Attachment:</p> <p>Notes:</p> |     |    |
| P-12                     | <p>The time it takes the inside of the cabinet to achieve the setpoint temperature measured in the absolute centre, must not exceed the time-frame specified in the appendix.</p> <p>Have the criteria been met?</p> <p>Attachment:</p> <p>Notes:</p>   |     |    |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_  
 Conducted by: \_\_\_\_\_  
 Inspected/verified by: \_\_\_\_\_  
 Model: \_\_\_\_\_ SN: \_\_\_\_\_

Deviations from the specifications dictated in the PQ, are to be reported in the deviation report.  
The PQ is concluded if all criteria of acceptance are approved and the possible deviations are rectified or accepted.

| Measurements – Hold-over |   |     |    |
|--------------------------|---|-----|----|
|                          |   | Yes | No |
| P-13                     | <p>The test is intended to provide substantiation for the time it takes for the temperature inside the cabinet to reach the end temperature specified in the appendix.<br/>Ambient temperature and setpoint temperature is specified in P-5.</p> <p>The temperature inside the cabinet must be stabilised – Where all the points in the working space have reached and maintained the same temperature throughout, the tolerances are specified in P-6.</p> <p>When the system is stable, turn off the power to the cabinet.</p> <p>The measurements, throughout the hold-over test, must be documented and attached the PQ.</p> <p>Attachment:</p> <p>Notes:</p> |     |    |
| P-14                     | <p>The times it takes the inside of the cabinet to reach the end temperature, must at least be the time specified in the appendix.</p> <p>Duration: _____</p> <p>Have the criteria been met?</p> <p>Attachment:</p> <p>Notes:</p>   |     |    |

Conducted by: \_\_\_\_\_ Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_

Inspected/verified by: \_\_\_\_\_

Model: \_\_\_\_\_ SN: \_\_\_\_\_

---

## Deviation Report

Deviations to the criteria of acceptance are to be documented in the deviation report. A separate deviation report shall be made for each deviation. Mark the entry with the relevant "P-ID" specified in the left column in the test specifications.

P-ID: \_\_\_\_\_

### Description of deviation:

---

---

---

---

### Extent to which the deviation has been alleviated:

---

---

---

---

---

### Additional notes:

---

---

---

---

---

Person responsible for test:

Name: \_\_\_\_\_

Date: \_\_\_\_\_

Organisation: \_\_\_\_\_

Signature: \_\_\_\_\_

Person responsible for verification of test:

Name: \_\_\_\_\_

Date: \_\_\_\_\_

Organisation: \_\_\_\_\_

Signature: \_\_\_\_\_

Model: \_\_\_\_\_

SN: \_\_\_\_\_

Approval of the test results – Performance Qualification (PQ)

- ☐ The steps in the Performance Qualification – PQ were completed with positive results
- ☐ The steps in the Performance Qualification – PQ were completed with negative results

ID of steps with negative results: \_\_\_\_\_

---

---

---

---

---

Organisation/Responsible party: \_\_\_\_\_

Trainer/Responsible party: \_\_\_\_\_

Stamp & Signature

Stamp & Signature

Tel.

Tel.

e-mail

e-mail

Location & Date

Location & Date

Model: \_\_\_\_\_ SN: \_\_\_\_\_

## Appendix

|   | Tolerances | Door opening recovery time (minutes) *** | Permissible deviation | Pull-down (minutes) | Permissible deviation | Hold-over range* | Hold-over | Permissible deviation |
|---|------------|--|-----------------------|---------------------|-----------------------|------------------|-----------|-----------------------|
| <b>Refrigerators (RR)</b>                             |            |  |                       |                     |                       |                  |           |                       |
| <b>BioMidi</b>  |            |  |                       |                     |                       |                  |           |                       |
| 425 (Solid door)                                      | +/- 3K     | 3 min.                                   | -                     | 20 min.             | -                     | -5 °C → 10 °C    | 63 min.   | -                     |
| 425 (Glass door)                                      | +/- 3K     | 4 min.                                   | -                     | 25 min.             | -                     | -5 °C → 10 °C    | 37 min.   | -                     |
| 625 (Solid door)                                      | +/- 3K     | 3 min.                                   | -                     | 20 min.             | -                     | -5 °C → 10 °C    | 63 min.   | -                     |
| 625 (Glass door)                                      | +/- 3K     | 4 min.                                   | -                     | 25 min.             | -                     | -5 °C → 10 °C    | 37 min.   | -                     |
| <b>Refrigerators with extended refrigeration (ER)</b> |            |  |                       |                     |                       |                  |           |                       |
| <b>BioPlus</b>  |            |  |                       |                     |                       |                  |           |                       |
| 500 (Solid door)                                      | +/- 2K     | 3 min.                                   | -                     | 22 min.             | -                     | -5 °C → 10 °C    | 72 min.   | -                     |
| 500 (Glass door)                                      | +/- 2K     | 4 min.                                   | -                     | 28 min.             | -                     | -5 °C → 10 °C    | 42 min.   | -                     |
| 600D/600W (Solid door)                                | +/- 2K     | 3 min.                                   | -                     | 20 min.             | -                     | -5 °C → 10 °C    | 70 min.   | -                     |
| 600D/600W (Glass door)                                | +/- 2K     | 4 min.                                   | -                     | 25 min.             | -                     | -5 °C → 10 °C    | 41 min.   | -                     |
| 660D/660W (Solid door)                                | +/- 2K     | 3 min.                                   | -                     | 20 min.             | -                     | -5 °C → 10 °C    | 70 min.   | -                     |
| 660D/660W (Glass door)                                | +/- 2K     | 4 min.                                   | -                     | 25 min.             | -                     | -5 °C → 10 °C    | 41 min.   | -                     |
| 930 (Solid door)                                      | +/- 2K     | 5 min.                                   | -                     | 22 min.             | -                     | -5 °C → 10 °C    | 65 min.   | -                     |
| 1270/1400 (Solid door)                                | +/- 2K     | 5 min.                                   | -                     | 23 min.             | -                     | -5 °C → 10 °C    | 78 min.   | -                     |
| 1270/1400 (Glass door)                                | +/- 2K     | 7 min.                                   | -                     | 29 min.             | -                     | -5 °C → 10 °C    | 45 min.   | -                     |

\*): The temperature span between the initial temperature and the end temperature in the hold-over test P-13,14.

\*\*): Please contact your local distributor for current information.

\*\*\*): 90° opening 1 minute

|       | Ambient Temperature | Setpoint temperature |
|-------|---------------------|----------------------|
| RR/ER | +25 °C              | +5 °C                |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_

Conducted by: \_\_\_\_\_

Inspected/verified by: \_\_\_\_\_

Model: \_\_\_\_\_ SN: \_\_\_\_\_



## Appendix

|                               | Tolerances | Door opening recovery time (minutes) *** | Permissible deviation | Pull-down (minutes) | Permissible deviation | Hold-over range* | Hold-over | Permissible deviation |
|-------------------------------|------------|--|-----------------------|---------------------|-----------------------|------------------|-----------|-----------------------|
| <b>FREEZERS (RF)</b>          |            |  |                       |                     |                       |                  |           |                       |
| <b>BioPlus</b>                |            |  |                       |                     |                       |                  |           |                       |
| 500                           | +/- 5K     | 7 min.                                   | -                     | 45 min.             | -                     | -20 °C → 10 °C   | 55 min.   | -                     |
| 600D/600W                     | +/- 5K     | 7 min.                                   | -                     | 42 min.             | -                     | -20 °C → 10 °C   | 55 min.   | -                     |
| 660D/660W                     | +/- 5K     | 7 min.                                   | -                     | 42 min.             | -                     | -20 °C → 10 °C   | 55 min.   | -                     |
| 930                           | +/- 5K     | -  | -                     | 76 min.             | -                     | -20 °C → 10 °C   | 63 min.   | -                     |
| 1270/1400                     | +/- 5K     | 10 min.                                  | -                     | 45 min.             | -                     | -20 °C → 10 °C   | 58 min.   | -                     |
| <b>EXTENDED FREEZERS (EF)</b> |            |  |                       |                     |                       |                  |           |                       |
| <b>BioMidi</b>                |            |  |                       |                     |                       |                  |           |                       |
| 425                           | +/- 9K     | 40 min.                                  | -                     | 107 min.            | -                     | -40 °C → 10 °C   | 108 min.  | -                     |
| <b>BioPlus</b>                |            |  |                       |                     |                       |                  |           |                       |
| 600W/660W                     | +/- 10K    | 30 min.                                  | -                     | 215 min.            | -                     | -35 °C → 10 °C   | 170 min.  | -                     |

\*): The temperature span between the initial temperature and the end temperature in the hold-over test P-13,14.

\*\*): Please contact your local distributor for current information.

\*\*\*): 90° opening 1 minute

|                | Ambient Temperature | Setpoint temperature |
|----------------|---------------------|----------------------|
| RF             | +25 °C              | -20 °C               |
| EF (425)       | +25 °C              | -40 °C               |
| EF (600W/660W) | +25 °C              | -35 °C               |

Name: \_\_\_\_\_ Signature: \_\_\_\_\_ Approved (Yes/No): \_\_\_\_\_ Date: \_\_\_\_\_

Conducted by: \_\_\_\_\_

Inspected/verified by: \_\_\_\_\_

Model: \_\_\_\_\_ SN: \_\_\_\_\_

---



**Gram Scientific Aps**

Aage Grams Vej 1 · 6500 Vojens · Denmark

Tel: +45 73 20 13 00

e-mail: [info@gram-bioline.com](mailto:info@gram-bioline.com)

[www.gram-bioline.com](http://www.gram-bioline.com)



Biostorage you can depend on